

Endobon® Xenograft Granules With OsseoGuard® And OsseoGuard Flex® Barrier Membranes

Tissue Management Treatment Solutions



Bone Graft Substitute

- Bovine-derived hydroxyapatite that has been fully deproteinized by a two-step, high temperature process for safety.
- An essentially non-resorbable material that is ideally suited for regeneration of bone defects when effective space maintenance is required.
- Osseointegrative due to the interconnecting micro and macro pores for bony integration, which facilitate graft stability and vascular ingrowth¹.
- Packaged in easy to open dishes. Large volumes (5 ml and 8 ml) are individually packaged in 1 ml containers for sterility.



Endobon Xenograft Granules adhere to one another when hydrated for easy transfer to the defect.

Manufacturer: Biomet France Sarl

Small Granules





500–1000 μm particle size typically preferred for grafting smaller defects, such as in extraction sockets.

Large Granules



1000–2000 μm particle size typically preferred for grafting large defects, such as sinus elevations because less material is needed with larger-sized particles.



SEM images of Endobon Xenograft Granules at 20x and 100x showing the micro and macro pores in the particles.

Endobon Xenograft Granules Are Indicated For Dental And/Or Oral Surgical Procedures, Such As:

- Alveolar ridge augmentation/reconstruction
- Filling of bone defects after root resection, cystectomy and apicectomy
- Filling socket after tooth extraction
- Sinus elevation

Barrier Membranes

- Resorbable collagen membranes designed for optimal strength, resorption, handling and biocompatibility.
- Made of highly purified collagen from safe bovine sources.
- A proprietary manufacturing process provides both membranes with a long resorption profile (6–9 months); well suited for Guided Bone Regeneration (GBR) procedures.
- OsseoGuard Membranes provide a protective barrier against soft-tissue invasion of a defect space.
- Two different levels of flexibility for ease of use in various clinical scenarios.
- Ability to tack or suture if desired.
- Three sizes are supplied sterile and are in double peel pouches for different defect sizes.



Manufacturer: Collagen Matrix, Inc., Oakland, NJ



OsseoGuard

Slightly more rigid for space maintenance.



OsseoGuard Flex

Intact tissue membrane for a higher degree of flexibility.
Performs when primary closure has not been achieved.^{2*}




Indications:

- Extraction sockets
- Localized ridge augmentation
- Alveolar ridge reconstruction
- GBR in dehiscence defects
- GTR in periodontal defects



* Clinical experience with OsseoGuard Flex has shown that the membrane shows no signs of inflammation or infection in cases where primary closure has not been achieved. The exposed area is healed by soft tissue covering the exposure within a few weeks while the membrane maintains its barrier function.

2. Clinical cases on file with Zimmer Biomet Dental



Cross section of the OsseoGuard Membrane at 200x magnification.



Posterior Mandible Recent Extraction Defects




Fig. 1: Clinical appearance of the surgical site at the time of implant placement four weeks after tooth extraction.



Fig. 2 & 3: The osseous defects were grafted with autogenous bone and Endobon Xenograft Small Granules. The surgical site was covered with an OsseoGuard 20 x30 mm Resorbable Collagen Membrane.



- The OsseoGuard Membrane is designed for combined strength, resorption and handling.
- Made of highly purified Type I collagen, derived from bovine Achilles tendon.

This provides:

- Combined strength to support suturing and good handling characteristics.
- A long resorption profile (6–9 months) suited for the healing time required in many GBR procedures.³

Manufacturer: Collagen Matrix, Inc.,
Oakland, NJ



Fig. 4: The surgical site was closed with sutures.



Fig. 5: Clinical appearance one month post-implant placement. Good epithelialization of the soft tissue is observed.



Fig. 6: Three months post-implant placement, the soft tissue has healed completely. The implants are ready for second stage surgery and healing abutment connection.



Fig. 7: Placement of the definitive restoration five months post-surgery.



Fig. 8: Clinical appearance nine months post-surgery. Note the healthy soft tissues.



Fig. 9: Periapical radiograph taken nine months post-surgery. Note the regenerated bone and graft integration.

Clinical Images Provided By: Dr. Francisco Enrile, Huelva, Spain.

3. Yuen D, Ulrich JB, Zuclich G, Homg-Ban L, Li S. Prediction of in vivo stability of a resorbable, reconstituted type I collagen membrane by in vitro methods. Society for Biomaterials, 2000.




Image illustrates the strength of the OsseoGuard Flex Membrane.

Maxillary Molar Post-Extraction Defects

Clinical images originally published by Dr. Robert del Castillo.[†]



Fig. 1: Extraction socket of first maxillary molar.



Fig. 2: Extraction socket grafted with Endobon Xenograft Small Granules and covered with an OsseoGuard Flex Membrane.




Fig. 3: The edges of the membrane were positioned under the soft tissue and secured with resorbable sutures.



Fig. 4: Healing was uneventful. The soft tissue was epithelializing over the OsseoGuard Flex Membrane two weeks postoperatively.



Fig. 5: The site was completely covered four weeks after the extraction.



Fig. 6: At four months postoperatively, a radiograph of the graft site showed excellent containment of the graft material.



Fig. 7: At four months postoperatively, the socket was healed and ready for implant placement.



Fig. 8: A 6.0 mm diameter Zimmer Biomet Dental Implant with a 5.0 mm platform was placed four months postoperatively.



Fig. 9: The implant was left sub-merged for two months of healing.

- The OsseoGuard Flex Membrane performs when primary closure has not been achieved.^{2*}
- The OsseoGuard Flex Membrane is designed for combined strength and drapability, resorption and handling.
- Made from highly purified intact bovine dermis tissue composed of Type I and Type III collagen.

This provides:

- Flexibility to drape over the defects.
- A long resorption profile (6–9 months) suited for the healing time required in many GBR procedures.⁴
- The ability to aid in gingival healing even when left exposed in a posterior molar extraction site.^{5**}

Manufacturer: Collagen Matrix, Inc.,
Oakland, NJ

** Primary closure is recommended. If exposed, resorption time will be shorter.

4. Yuen D, Ulreich JB, Zuclich G, Homg-Ban L, Li S. Prediction of in vivo stability of a resorbable, reconstituted Type 1 collagen membrane by in vitro methods. Society for Biomaterials, 2000.

5. del Castillo R[†]. Grafting of an extracted maxillary first-molar socket using a new, flexible resorbable collagen membrane for ridge preservation in advance of implant placement. Inside Dentistry, October 2011. 94-96.

[†] Dr. del Castillo has a financial relationship with Zimmer Biomet Dental resulting from speaking engagements, consulting engagements and other retained services.

Anterior Ridge Augmentation




Fig. 1: Extraction sockets of the four maxillary incisors and immediate implant placement.



Fig. 2: Grafting with Endobon Xenograft Small Granules covered by an OsseoGuard Resorbable Collagen Membrane.



Fig. 3: The soft-tissue flaps were closed and sutured.



Fig. 4: Clinical appearance of soft tissue showing excellent soft-tissue healing after four months.




Fig. 5: Clinical appearance of the regenerated site at four months after removing the remnants of the membrane.



Fig. 6: Occlusal view after four months.

Post-Extraction Defects in the Aesthetic Zone



Fig. 7: Post-extraction defects in the maxilla right central and lateral incisor area.



Fig. 8: Occlusal view of the extraction site defects.



Fig. 9: Facial view of dehiscence defects after implant placement.



Fig. 10: Occlusal view of implants and defects.



Fig. 11: Grafting with Endobon Xenograft Small Granules covered by an OsseoGuard Resorbable Collagen Membrane.



Fig. 12: Regeneration at four months after removing the remnants of the membrane.

Clinical Images Provided By: Dr. Xavier Vela[†], Barcelona, Spain

[†] Dr. Xavier Vela has a financial relationship with Zimmer Biomet Dental resulting from speaking engagements, consulting engagements and other retained services.

Histological Study of Endobon Xenograft Granules in Sinus Floor Augmentation



“The Clinical and Histological Efficacy of Xenograft Granules for Maxillary Sinus Floor Augmentation”

A study led by Dr. Myron Nevins[†] at the Harvard School of Dentistry published in The International Journal Of Periodontics & Restorative Dentistry (2011 Jun;31(3):227-235), highlights the positive results that clinicians achieved when using Endobon Xenograft Granules in patients requiring sinus augmentation procedures prior to implant placement.

At six months postoperatively, the following observations were made:

- Bone formation at the osteotomy site ranging from 16.2% to 43.6% was observed in all patients.
- Histologic evaluation showed Endobon Xenograft Granules to be integrated and surrounded by woven bone and in close contact with the particles.
- No inflammatory cells were present and there were no signs of Xenograft resorption.
- Evidence was observed of woven bone undergoing remodeling and maturing to well-organized lamellar bone.

Some areas of the newly formed bone were undergoing remodeling, maturing from woven bone (WB) to well-organized lamellar bone (LB). XG: Xenograft Granules.



Scientific References

- Barone A, Todisco M, Ludovichetti M, Gualini F, Torres-Lagares D, Aggstaller H. A prospective randomized-controlled study of Endobon used in extraction sites: a clinical and histological evaluation. *Int J Periodontics Restorative Dent* 2013;33:795–802. doi: 10.11607/prd.1690. (Endobon and OsseoGuard)
- Ramírez-Fernández MP, Calvo-Guirado JL, Delgado-Ruiz RA, Maté-Sánchez Del Val JE, Negri B, Peñarrocha Diago M. Ultrastructural study by back scattered electron imaging and elemental microanalysis of biomaterial-to-bone interface and mineral degradation of bovine xenografts in maxillary sinus floor elevation. *Clin Oral Implants Res*. 2013;24:645–651. doi: 10.1111/j.1600-0501.2012.02439. (Endobon and OsseoGuard)
- Testori T[†], Izzi G, Manzoni L, Fratto G, Piatelli A, Weinstein RL. High temperature-treated bovine porous hydroxyapatite in sinus augmentation procedures: a case report. *Int J Periodontics Restorative Dent*. 2012; 32(3):295-301. (Endobon)
- Mazor Z. A prospective, randomized-controlled study of implants placed simultaneously in maxillary sinus augmentation cases: a clinical and histological evaluation. 27th Annual Meeting Academy of Osseointegration. 2012, March 1-3:Phoenix, Arizona, USA. www.iird.com/pdf/p14-mazor.pdf. (Endobon and OsseoGuard)
- Weinländer M[†], Krennmair G, Schmidinger S, Plenck H, Piatelli A. The Mirror Study. A prospective randomized controlled study of Endobon used in maxillary sinus augmentation. Histological evaluation. 20th Annual Scientific Meeting of the European Association of Osseointegration. 2011, October 13-15: Athens, Greece. Abstract 341 *Clinical Oral Implants Res* 2011;22(9): p 1033. (supplement) (Endobon and OsseoGuard)
- De Angelis N[†], Felice P, Pellegrino G, Camurati A, Gambino P, Esposito M[†]. Guided bone regeneration with and without a bone substitute at single post-extractive implants: 1-year post-loading results from a pragmatic multicenter randomised controlled trial. *Eur J Oral Implantol*. 2011;4(4):313-25. (Endobon and OsseoGuard)
- Nevins M, Camelo M, De Angelis N, Hanratty JJ, Khang WG, Kwon J-J[†], Rasperini G, Rocchietta D, Schupbach P, Kim DM[†]. The clinical and histologic efficacy of xenograft granules for maxillary sinus floor augmentation. *Int J Periodontics Restorative Dent*. 2011;31:227-235. (Endobon and OsseoGuard)
- Ramírez-Fernández MR, Guirado JL, Ruiz RA, Sánchez de-Val JE, Ortega VV, Olmos LM. Bone response to hydroxyapatites with open porosity of animal origin (porcine [OsteoBio[®] mp3] and bovine [Endobon]): a radiological and histomorphometric study. *Clin Oral Implants Res*. 2011;22(7):767-73. (Endobon)
- Ramírez-Fernández MaP, Calvo-Guirado JL, Arcesio-Delgado Ruiz R, Maté-Sánchez de Val JE, Gómez-Moreno G, Guardia J. Experimental model of bone response to xenografts of bovine origin (Endobon): a radiological and histomorphometric study. *Clin Oral Implants Res* 2011;22(7):727-34. (Endobon)
- del Castillo R[†]. Grafting of an Extracted Maxillary First-Molar Socket using a New, Flexible Resorbable Collagen Membrane for Ridge Preservation in Advance of Implant Placement. *Inside Dentistry*, October 2011;94-96. www.dentalaegis.com/id. (Endobon and OsseoGuard Flex)
- Spies CK, Schnürer S, Gotterbarm T, Breusch SJ. Efficacy of Bone Source[™] and Cementek[™] in comparison with Endobon in critical size metaphyseal defects, using a minipig model. *J Appl Biomater Biomech*. 2010 Sep-Dec;8(3):175-85. (Endobon)
- Aguilar-Salvatierra A, Gómez-Moreno G, Martín-Piedra M.Á., M. Cabrera M, López-Gallardo C, Guardia J, Castillo T, Calvo-Guirado JL, Ramirez-Fernandez MP, Delgado-Ruiz R, And López-Mari L. Radiological and histomorphometric analysis of two xenografts: an experimental study. International Association of Dental Research/ General Session 2010. July 14-17, Barcelona, Spain. <https://iadr.confex.com/iadr/2010barce/webprogram/Session23668.html>. (Endobon)
- Schnettler R, Knöss PD, Heiss C, Stahl JP, Meyer C, Kilian O, Wenisch S, Alt V. Enhancement of bone formation in hydroxyapatite implants by rhBMP-2 coating. *J Biomed Mater Res B Appl Biomater*. 2009;90(1):75-81. (Endobon)
- Jensen SS, Aaboe M, Pinholt EM, Hjørtting-Hansen E, Melsen F, Ruyter IE. Tissue reaction and material characteristics of four bone substitutes. *Int J Oral Maxillofac Implants*. 1996;11(1):55-66. (Endobon)
- Wiltfang J, Merten HA, Wiltfang J. Ectopic bone formation with the help of growth factor bFGF. *J Craniomaxillofac Surg*. 1996 Oct;24(5):300-4. (Endobon)

[†] These clinicians have current or past financial relationships with Zimmer Biomet Dental resulting from speaking engagements, consulting engagements and other retained services.

Ordering Information

OsseoGuard And OsseoGuard Flex Barrier Membranes

| Size (mm) | OsseoGuard Membrane | OsseoGuard Flex Membrane |
|------------|---------------------|--------------------------|
| 15 x 20 mm | OG1520 | OGF1520 |
| 20 x 30 mm | OG2030 | OGF2030 |
| 30 x 40 mm | OG3040 | OGF3040 |

Shelf Life: 3 Years



Manufacturer: Collagen Matrix, Inc., Oakland, NJ

Endobon Xenograft Granules

| Volume (ml) | Small Granules 500–1,000 µm | Large Granules 1,000–2,000 µm |
|-------------|--------------------------------|----------------------------------|
| 0.5 ml | ROX05 | N/A |
| 1 ml | ROX10 | N/A |
| 2 ml | ROX20 | ROXLG20 |
| 5 ml | N/A | ROXLG50 |
| 8 ml | N/A | ROXLG80 |

Shelf-Life: 18 Months



Manufacturer: Biomet France Sarl.

**Want To Facilitate Re-Growth Of Bone And Soft Tissue?
Consider Zimmer Biomet Dental Regenerative Options.**

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Tel: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). OsseoGuard and OsseoGuard Flex Membranes are manufactured by Collagen Matrix, Inc. Collagen Matrix, Inc. is not a subsidiary of Zimmer Biomet Holdings, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0042 REV A 03/18 ©2018 Zimmer Biomet. All rights reserved.

