

# Certain® Implant System

## Product Catalog



#### **Important Product Information**

##### **Biomet 3i Dental Implants**

For detailed product information on all Biomet 3i Dental Implants, refer to the Biomet 3i Dental Implant IFU (P-IIS086GI) located at [ifu.biomet3i.com](http://ifu.biomet3i.com).

##### **Biomet 3i Restorative Products**

For detailed product information on all Biomet 3i Restorative Products, refer to the Biomet 3i Restorative Products IFU (P-IIS086GR) located at [ifu.biomet3i.com](http://ifu.biomet3i.com).

# Table Of Contents

## Implants

### Tapered Certain Internal Connection Implants

|                                          |   |
|------------------------------------------|---|
| T3® With DCD® Platform Switched          | 4 |
| T3 With DCD Non-Platform Switched        | 4 |
| T3 Platform Switched                     | 4 |
| T3 Non-Platform Switched                 | 4 |
| Full Osseotite® Tapered Certain PREVAIL® | 4 |
| Full Osseotite Tapered Certain           | 4 |

### Parallel Walled Certain Internal Connection Implants

|                                     |   |
|-------------------------------------|---|
| T3 With DCD Platform Switched       | 5 |
| T3 With DCD Non-Platform Switched   | 5 |
| T3 Platform Switched                | 5 |
| T3 Non-Platform Switched            | 5 |
| Osseotite 2 Parallel Walled Certain | 5 |

## Site Preparation Technologies

### Tapered Navigator® System For Guided Surgery

|                                  |   |
|----------------------------------|---|
| Tapered Navigator Surgical Kit   | 6 |
| Tapered Navigator Tap Kit        | 8 |
| Tapered Navigator Laboratory Kit | 9 |

### Parallel Walled Navigator System For Guided Surgery

|                                          |    |
|------------------------------------------|----|
| Parallel Walled Navigator Surgical Kit   | 10 |
| Parallel Walled Navigator Laboratory Kit | 12 |

### Surgical Kits

|                               |    |
|-------------------------------|----|
| Tapered Implant Surgical Kits | 13 |
| Parallel Walled Surgical Kits | 16 |

### Surgical Drills

|                                                  |    |
|--------------------------------------------------|----|
| W & H Surgical Drilling Units                    | 18 |
| Advanced Cutting Technology® (ACT®) Twist Drills | 19 |
| Quad Shaping Drills For Tapered Implants         | 19 |
| Miscellaneous Drills                             | 19 |

### Surgical Components

|                                              |    |
|----------------------------------------------|----|
| Dense Bone Taps For Parallel Walled Implants | 20 |
| Dense Bone Taps For Tapered Implants         | 20 |
| Tapered Implant Depth Indicators             | 21 |
| Direction Indicators                         | 21 |
| Bone Profilers                               | 22 |
| Tapered Osteotomes                           | 22 |
| Summer Osteotomes                            | 22 |
| Sinus Elevation                              | 23 |

|               |    |
|---------------|----|
| Trephine Burs | 23 |
|---------------|----|

|                               |    |
|-------------------------------|----|
| Miscellaneous Instrumentation | 24 |
|-------------------------------|----|

## Restorative Components

### Healing Abutments

|                   |    |
|-------------------|----|
| BellaTek® Encode® | 27 |
| EP® One-piece     | 27 |

### Provisional Restorations

|                     |    |
|---------------------|----|
| PreFormance® Posts  | 28 |
| Temporary Cylinders | 28 |

### Definitive Restorations

|                                            |    |
|--------------------------------------------|----|
| GingiHue® Posts                            | 29 |
| UCLA Abutments                             | 29 |
| GingiHue Post And UCLA Abutment Components | 30 |
| Low Profile Abutments                      | 31 |
| Low Profile Abutment Components            | 32 |

### Overdenture Restorations

|                         |    |
|-------------------------|----|
| LOCATOR R-Tx® Abutments | 33 |
| LOCATOR® Abutments      | 33 |

### Screws

|                                 |    |
|---------------------------------|----|
| Abutment Screws                 | 34 |
| Retaining Screws                | 34 |
| Waxing Screws And Try-In Screws | 34 |

### BellaTek Bars And Frameworks

### Mechanical Driver Kits

### Prosthodontic Screwdriver

### Drivers, Driver Tips And Miscellaneous Restorative Products


# Tapered Certain Internal Connection

4

## Implants

### T3 With DCD Platform Switched Implants


Commercially Pure Titanium



| Length  | 4/3 mmP  | 5/4 mmP  | 6/5 mmP  |
|---------|----------|----------|----------|
| 8.5 mm  | BNPT4385 | BNPT5485 | BNPT6585 |
| 10 mm   | BNPT4310 | BNPT5410 | BNPT6510 |
| 11.5 mm | BNPT4311 | BNPT5411 | BNPT6511 |
| 13 mm   | BNPT4313 | BNPT5413 | BNPT6513 |
| 15 mm   | BNPT4315 | BNPT5415 | BNPT6515 |

### T3 With DCD Non-Platform Switched Implants

Commercially Pure Titanium



| 3.25 mmD | 4 mmD   | 5 mmD   | 6 mmD   |
|----------|---------|---------|---------|
| BNST3285 | BNST485 | BNST585 | BNST685 |
| BNST3210 | BNST410 | BNST510 | BNST610 |
| BNST3211 | BNST411 | BNST511 | BNST611 |
| BNST3213 | BNST413 | BNST513 | BNST613 |
| BNST3215 | BNST415 | BNST515 | BNST615 |

### T3 Platform Switched Implants

Commercially Pure Titanium

| Length  | 4/3 mmP  | 5/4 mmP  | 6/5 mmP  |
|---------|----------|----------|----------|
| 8.5 mm  | BOPT4385 | BOPT5485 | BOPT6585 |
| 10 mm   | BOPT4310 | BOPT5410 | BOPT6510 |
| 11.5 mm | BOPT4311 | BOPT5411 | BOPT6511 |
| 13 mm   | BOPT4313 | BOPT5413 | BOPT6513 |
| 15 mm   | BOPT4315 | BOPT5415 | BOPT6515 |


### T3 Non-Platform Switched Implants

Commercially Pure Titanium

| 3.25 mmD | 4 mmD   | 5 mmD   | 6 mmD   |
|----------|---------|---------|---------|
| BOST3285 | BOST485 | BOST585 | BOST685 |
| BOST3210 | BOST410 | BOST510 | BOST610 |
| BOST3211 | BOST411 | BOST511 | BOST611 |
| BOST3213 | BOST413 | BOST513 | BOST613 |
| BOST3215 | BOST415 | BOST515 | BOST615 |

### Full Osseotite Tapered Certain PREVAIL


Commercially Pure Titanium



| Length  | 4/3 mmP   | 5/4 mmP   | 6/5 mmP   |
|---------|-----------|-----------|-----------|
| 8.5 mm  | XIITP4385 | XIITP5485 | XIITP6585 |
| 10 mm   | XIITP4310 | XIITP5410 | XIITP6510 |
| 11.5 mm | XIITP4311 | XIITP5411 | XIITP6511 |
| 13 mm   | XIITP4313 | XIITP5413 | XIITP6513 |
| 15 mm   | XIITP4315 | XIITP5415 | XIITP6515 |

### Full Osseotite Tapered Certain

Commercially Pure Titanium



| 3.25 mmD  | 4 mmD    | 5 mmD    | 6.0 mmD  |
|-----------|----------|----------|----------|
| XIFNT3285 | XIFNT485 | XIFNT585 | XIFNT685 |
| XIFNT3210 | XIFNT410 | XIFNT510 | XIFNT610 |
| XIFNT3211 | XIFNT411 | XIFNT511 | XIFNT611 |
| XIFNT3213 | XIFNT413 | XIFNT513 | XIFNT613 |
| XIFNT3215 | XIFNT415 | XIFNT515 | XIFNT615 |

| Cover Screw Flat (included) | IMCSF34 | ICSF41 | ICSF50 |
|-----------------------------|---------|--------|--------|
|                             |         |        |        |

| IMCSF34 | ICSF41 | ICSF50 | ICSF60 |
|---------|--------|--------|--------|
|         |        |        |        |

\*Non-Flared 4.1 mmP Cover Screw ICS400 is also available.

# Parallel Walled Certain Internal Connection

5


## Implants

### T3 With DCD Platform Switched Implants

Commercially Pure Titanium




Contemporary Hybrid Surface



| Length  | 4/3 mmP  | 5/4 mmP  | 6/5 mmP  |
|---------|----------|----------|----------|
| 8.5 mm  | BNPS4385 | BNPS5485 | BNPS6585 |
| 10 mm   | BNPS4310 | BNPS5410 | BNPS6510 |
| 11.5 mm | BNPS4311 | BNPS5411 | BNPS6511 |
| 13 mm   | BNPS4313 | BNPS5413 | BNPS6513 |
| 15 mm   | BNPS4315 | BNPS5415 | BNPS6515 |
| 18 mm   | —        | —        | —        |

### T3 With DCD Non-Platform Switched Implants

Commercially Pure Titanium



| Length  | 3.25 mmD | 4 mmD   | 5 mmD   | 6 mmD   |
|---------|----------|---------|---------|---------|
| 8.5 mm  | BNSS385  | BNSS485 | BNSS585 | BNSS685 |
| 10 mm   | BNSS310  | BNSS410 | BNSS510 | BNSS610 |
| 11.5 mm | BNSS311  | BNSS411 | BNSS511 | BNSS611 |
| 13 mm   | BNSS313  | BNSS413 | BNSS513 | BNSS613 |
| 15 mm   | BNSS315  | BNSS415 | BNSS515 | BNSS615 |
| 18 mm   | BNSS318  | BNSS418 | —       | —       |

### T3 Platform Switched Implants

Commercially Pure Titanium

| Length  | 4/3 mmP  | 5/4 mmP  | 6/5 mmP  |
|---------|----------|----------|----------|
| 8.5 mm  | BOPS4385 | BOPS5485 | BOPS6585 |
| 10 mm   | BOPS4310 | BOPS5410 | BOPS6510 |
| 11.5 mm | BOPS4311 | BOPS5411 | BOPS6511 |
| 13 mm   | BOPS4313 | BOPS5413 | BOPS6513 |
| 15 mm   | BOPS4315 | BOPS5415 | BOPS6515 |
| 18 mm   | —        | —        | —        |

### T3 Non-Platform Switched Implants

Commercially Pure Titanium


| Length  | 3.25 mmD | 4 mmD   | 5 mmD   | 6 mmD   |
|---------|----------|---------|---------|---------|
| 8.5 mm  | BOSS385  | BOSS485 | BOSS585 | BOSS685 |
| 10 mm   | BOSS310  | BOSS410 | BOSS510 | BOSS610 |
| 11.5 mm | BOSS311  | BOSS411 | BOSS511 | BOSS611 |
| 13 mm   | BOSS313  | BOSS413 | BOSS513 | BOSS613 |
| 15 mm   | BOSS315  | BOSS415 | BOSS515 | BOSS615 |
| 18 mm   | BOSS318  | BOSS418 | —       | —       |

### Osseotite 2 Parallel Walled Certain

Commercially Pure Titanium



Full Osseotite Surface



| Length  | 3.25 mmD  | 4 mmD     | 5 mmD     | 6 mmD     |
|---------|-----------|-----------|-----------|-----------|
| 8.5 mm  | XIFOSM385 | XIFOSS485 | XIFOSS585 | XIFOSS685 |
| 10 mm   | XIFOSM310 | XIFOSS410 | XIFOSS510 | XIFOSS610 |
| 11.5 mm | XIFOSM311 | XIFOSS411 | XIFOSS511 | XIFOSS611 |
| 13 mm   | XIFOSM313 | XIFOSS413 | XIFOSS513 | XIFOSS613 |
| 15 mm   | XIFOSM315 | XIFOSS415 | XIFOSS515 | XIFOSS615 |

| Cover Screw Flat (included) | IMCSF34 | ICSF41 | ICSF50 | ICSF60 |
|-----------------------------|---------|--------|--------|--------|
|                             |         |        |        |        |

\*Non-Flared 4.1 mmP Cover Screw ICS400 is also available.

# Tapered Navigator System For Guided Surgery

Site Preparation Technologies

## Tapered Navigator Surgical Kit

### For Tapered Implants Including:

- Certain 3.25, 4 and 5 mm Implants
- Prevail 4/3 and 5/4 Implants


The Navigator System helps clinicians transform computer-based planning into practical implant placement with an opportunity for immediate provisionalization.

Features specific to the Navigator System include:

- Depth-specific instrumentation for preparation and implant placement
- Control of hex-orientation to fabricate and deliver provisional restorations immediately following implant placement
- Open architecture compatibility that provides freedom to choose among leading planning software solutions

SGTIKIT

Tapered Navigator Surgical Kit



| Description                                                  | Item No. | SGTIKIT |
|--------------------------------------------------------------|----------|---------|
| Tapered Navigator Surgical Tray                              | SGTTRAY  | •       |
| Tapered Navigator Tissue Punch, 4.1 mmD                      | SGTTP4   | •       |
| Tapered Navigator Tissue Punch, 5 mmD                        | SGTTP5   | •       |
| Tapered Navigator Cortical Perforator, 4.1 mmD - Short       | SGTCP4S  | •       |
| Tapered Navigator Cortical Perforator, 4.1 mmD - Long        | SGTCP4L  | •       |
| Tapered Navigator Cortical Perforator, 5 mmD - Short         | SGTCP5S  | •       |
| Tapered Navigator Cortical Perforator, 5 mmD - Long          | SGTCP5L  | •       |
| Tapered Navigator Twist Drill, 1.9 mmD x A(L)                | SGTD19A  | •       |
| Tapered Navigator Twist Drill, 1.9 mmD x B(L)                | SGTD19B  | •       |
| Tapered Navigator Twist Drill, 1.9 mmD x C(L)                | SGTD19C  | •       |
| Tapered Navigator Twist Drill, 1.9 mmD x D(L)                | SGTD19D  | •       |
| Tapered Navigator Twist Drill, 1.9 mmD x E(L)                | SGTD19E  | •       |
| Tapered Navigator Twist Drill, 1.9 mmD x F(L)                | SGTD19F  | •       |
| Tapered Navigator Twist Drill, 1.9 mmD x G(L)                | SGTD19G  | •       |
| Tapered Navigator Countersink Drill, 3.4 mmD - Short         | MSGTCD5  | •       |
| Tapered Navigator Countersink Drill, 4.1 mmD - Short         | SGTCD4S  | •       |
| Tapered Navigator Countersink Drill, 5 mmD - Short           | SGTCD5S  | •       |
| Tapered Navigator Countersink Drill, 3.4 mmD - Long          | MSGTCDL  | •       |
| Tapered Navigator Countersink Drill, 4.1 mmD - Long          | SGTCD4L  | •       |
| Tapered Navigator Countersink Drill, 5 mmD - Long            | SGTCD5L  | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 8.5 mmL - Short  | MSGTD85S | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 10 mmL - Short   | MSGTD10S | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 11.5 mmL - Short | MSGTD11S | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 13 mmL - Short   | MSGTD13S | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 15 mmL - Short   | MSGTD15S | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 8.5 mmL - Short   | SGTD485S | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 10 mmL - Short    | SGTD410S | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 11.5 mmL - Short  | SGTD411S | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 13 mmL - Short    | SGTD413S | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 15 mmL - Short    | SGTD415S | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 8.5 mmL - Short     | SGTD585S | •       |

## Tapered Navigator Surgical Kit (Continued)

| Description                                                 | Item No. | SGTIKIT |
|-------------------------------------------------------------|----------|---------|
| Tapered Navigator Shaping Drill, 5 mmD x 10 mmL - Short     | SGTD510S | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 11.5 mmL - Short   | SGTD511S | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 13 mmL - Short     | SGTD513S | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 15 mmL - Short     | SGTD515S | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 8.5 mmL - Long  | MSGTD85L | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 10 mmL - Long   | MSGTD10L | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 11.5 mmL - Long | MSGTD11L | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 13 mmL - Long   | MSGTD13L | •       |
| Tapered Navigator Shaping Drill, 3.25 mmD x 15 mmL - Long   | MSGTD15L | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 8.5 mmL - Long   | SGTD485L | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 10 mmL - Long    | SGTD410L | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 11.5 mmL - Long  | SGTD411L | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 13 mmL - Long    | SGTD413L | •       |
| Tapered Navigator Shaping Drill, 4.1 mmD x 15 mmL - Long    | SGTD415L | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 8.5 mmL - Long     | SGTD585L | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 10 mmL - Long      | SGTD510L | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 11.5 mmL - Long    | SGTD511L | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 13 mmL - Long      | SGTD513L | •       |
| Tapered Navigator Shaping Drill, 5 mmD x 15 mmL - Long      | SGTD515L | •       |
| Tapered Navigator Twist Drill Positioning Handle            | SGTDPHT  | •       |
| Tapered Navigator Reduction Drill Positioning Handle        | SGTDPHR  | •       |
| Handpiece Connector                                         | MDR10    | •       |
| Narrow Right Angle Large Hexed Driver Tip, 24 mmL           | RASH3N   | •       |
| Posterior Large Hexed Driver, 17 mmL                        | PHD02N   | •       |
| Ratchet Extension, 6 mmL - Short                            | RE100    | •       |
| Ratchet Extension, 15 mmL - Long                            | RE200    | •       |
| Ratchet Wrench                                              | WR150    | •       |
| Open End Wrench                                             | CW100    | •       |
| Navigator Certain Implant Mount, 3.4 mmD - Short            | MSGIIMS  | •       |
| Navigator Certain Implant Mount, 3.4 mmD - Long             | MSGIIML  | •       |
| Navigator Certain Implant Mount, 4.1 mmD - Short            | SGIIM4S  | •       |
| Navigator Certain Implant Mount, 4.1 mmD - Long             | SGIIM4L  | •       |
| Navigator Certain Implant Mount, 5 mmD - Short              | SGIIM5S  | •       |
| Navigator Certain Implant Mount, 5 mmD - Long               | SGIIM5L  | •       |
| Certain Manual Bone Profiler, 3.4 mmD x 4.1 mmP             | IMMBP    | •       |
| Certain Manual Bone Profiler, 4.1 mmD x 5 mmP               | IMBP4    | •       |
| Certain Manual Bone Profiler, 5 mmD x 6 mmP                 | IMBP5    | •       |

# Tapered Navigator System For Guided Surgery

Site Preparation Technologies

## Tapered Navigator Tap Kit

### For Tapered Implants Including:

- Certain 3.25, 4 and 5 mm Implants
- Prevail 4/3 and 5/4 Implants

### Tapered Navigator Tap Kit

(Included under SGTIKIT)



| Description                                             | Item No. | SGTIKIT |
|---------------------------------------------------------|----------|---------|
| Tapered Navigator Bone Tap Tray                         | SGTTAPT  | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 8.5 mmL - Short  | MSGTT85S | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 10 mmL - Short   | MSGTT10S | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 11.5 mmL - Short | MSGTT11S | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 13 mmL - Short   | MSGTT13S | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 15 mmL - Short   | MSGTT15S | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 8.5 mmL - Short   | SGTT485S | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 10 mmL - Short    | SGTT410S | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 11.5 mmL - Short  | SGTT411S | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 13 mmL - Short    | SGTT413S | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 15 mmL - Short    | SGTT415S | •       |
| Tapered Navigator Bone Tap, 5 mmD x 8.5 mmL - Short     | SGTT585S | •       |
| Tapered Navigator Bone Tap, 5 mmD x 10 mmL - Short      | SGTT510S | •       |
| Tapered Navigator Bone Tap, 5 mmD x 11.5 mmL - Short    | SGTT511S | •       |
| Tapered Navigator Bone Tap, 5 mmD x 13 mmL - Short      | SGTT513S | •       |
| Tapered Navigator Bone Tap, 5 mmD x 15 mmL - Short      | SGTT515S | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 8.5 mmL - Long   | MSGTT85L | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 10 mmL - Long    | MSGTT10L | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 11.5 mmL - Long  | MSGTT11L | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 13 mmL - Long    | MSGTT13L | •       |
| Tapered Navigator Bone Tap, 3.25 mmD x 15 mmL - Long    | MSGTT15L | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 8.5 mmL - Long    | SGTT485L | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 10 mmL - Long     | SGTT410L | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 11.5 mmL - Long   | SGTT411L | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 13 mmL - Long     | SGTT413L | •       |
| Tapered Navigator Bone Tap, 4.1 mmD x 15 mmL - Long     | SGTT415L | •       |
| Tapered Navigator Bone Tap, 5 mmD x 8.5 mmL - Long      | SGTT585L | •       |
| Tapered Navigator Bone Tap, 5 mmD x 10 mmL - Long       | SGTT510L | •       |
| Tapered Navigator Bone Tap, 5 mmD x 11.5 mmL - Long     | SGTT511L | •       |
| Tapered Navigator Bone Tap, 5 mmD x 13 mmL - Long       | SGTT513L | •       |
| Tapered Navigator Bone Tap, 5 mmD x 15 mmL - Long       | SGTT515L | •       |

## Tapered Navigator Laboratory Kit

### For Tapered Implants Including:

- Certain 3.25, 4 and 5 mm Implants
- Prevail 4/3 and 5/4 Implants

SGTILKIT

Tapered Navigator Laboratory Kit



| Description                                     | Item No. | SGTILKIT |
|-------------------------------------------------|----------|----------|
| Tapered Navigator Laboratory Tray               | SGTLTRAY | •        |
| Navigator Certain Analog Mount, 3.4 mmD - Short | MSGIAMS  | •        |
| Navigator Certain Analog Mount, 4.1 mmD - Short | SGIAM4S  | •        |
| Navigator Certain Analog Mount, 5 mmD - Short   | SGIAM5S  | •        |
| Navigator Certain Analog Mount, 3.4 mmD - Long  | MSGIAML  | •        |
| Navigator Certain Analog Mount, 4.1 mmD - Long  | SGIAM4L  | •        |
| Navigator Certain Analog Mount, 5 mmD - Long    | SGIAM5L  | •        |

# Parallel Walled Navigator System For Guided Surgery

Site Preparation Technologies

## Parallel Walled Navigator Surgical Kit

### For Parallel Walled Implants Including:

- Certain 3.25, 4 and 5 mm Implants

The Navigator System helps clinicians transform computer-based planning into practical implant placement with an opportunity for immediate provisionalization.

Features specific to the Navigator System include:

- Depth-specific instrumentation for preparation and implant placement
- Control of hex-orientation to fabricate and deliver provisional restorations immediately following implant placement
- Open architecture compatibility that provides freedom to choose among leading planning software solutions

SGKIT Parallel Walled Navigator Surgical Kit



| Description                            | Item No. | SGKIT |
|----------------------------------------|----------|-------|
| Navigator Surgical Tray                | SGTRAY   | •     |
| Navigator Tissue Punch, 4.1 mmD        | SGTP4    | •     |
| Navigator Tissue Punch, 5 mmD          | SGTP5    | •     |
| Navigator Starter Drill, 3.4 mmD       | MSGSD    | •     |
| Navigator Starter Drill, 3/4 mmD       | SGSD34   | •     |
| Navigator Starter Drill, 4.1 mmD       | SGSD4    | •     |
| Navigator Starter Drill, 4/5 mmD       | SGSD45   | •     |
| Navigator Certain Starter Drill, 5 mmD | SGISD5   | •     |
| Navigator Twist Drill, 2 mmD x A(L)    | SGD2A    | •     |
| Navigator Twist Drill, 2 mmD x B(L)    | SGD2B    | •     |
| Navigator Twist Drill, 2 mmD x C(L)    | SGD2C    | •     |
| Navigator Twist Drill, 2 mmD x D(L)    | SGD2D    | •     |
| Navigator Twist Drill, 2 mmD x E(L)    | SGD2E    | •     |
| Navigator Twist Drill, 2.75 mmD x A(L) | SGD275A  | •     |
| Navigator Twist Drill, 2.75 mmD x B(L) | SGD275B  | •     |
| Navigator Twist Drill, 2.75 mmD x C(L) | SGD275C  | •     |
| Navigator Twist Drill, 2.75 mmD x D(L) | SGD275D  | •     |
| Navigator Twist Drill, 2.75 mmD x E(L) | SGD275E  | •     |
| Navigator Twist Drill, 3 mmD x A(L)    | SGD3A    | •     |
| Navigator Twist Drill, 3 mmD x B(L)    | SGD3B    | •     |
| Navigator Twist Drill, 3 mmD x C(L)    | SGD3C    | •     |
| Navigator Twist Drill, 3 mmD x D(L)    | SGD3D    | •     |
| Navigator Twist Drill, 3 mmD x E(L)    | SGD3E    | •     |
| Navigator Twist Drill, 3.25 mmD x A(L) | SGD325A  | •     |
| Navigator Twist Drill, 3.25 mmD x B(L) | SGD325B  | •     |
| Navigator Twist Drill, 3.25 mmD x C(L) | SGD325C  | •     |
| Navigator Twist Drill, 3.25 mmD x D(L) | SGD325D  | •     |
| Navigator Twist Drill, 3.25 mmD x E(L) | SGD325E  | •     |
| Navigator Twist Drill, 3.85 mmD x A(L) | SGD385A  | •     |
| Navigator Twist Drill, 3.85 mmD x B(L) | SGD385B  | •     |
| Navigator Twist Drill, 3.85 mmD x C(L) | SGD385C  | •     |
| Navigator Twist Drill, 4.25 mmD x A(L) | SGD425A  | •     |

## Parallel Walled Navigator Surgical Kit (Continued)

| Description                                       | Item No. | SGKIT |
|---------------------------------------------------|----------|-------|
| Navigator Twist Drill, 4.25 mmD x B(L)            | SGD425B  | •     |
| Navigator Twist Drill, 4.25 mmD x C(L)            | SGD425C  | •     |
| Navigator Bone Tap, 3.25 mmD                      | MSGTAP   | •     |
| Navigator Bone Tap, 4.1 mmD                       | SGTAP4   | •     |
| Navigator Bone Tap, 4/5 mmD                       | SGTAP45  | •     |
| Navigator Bone Tap, 5 mmD                         | SGTAP5   | •     |
| Navigator Drill Positioning Handle 1              | SGDPH1   | •     |
| Navigator Drill Positioning Handle 2              | SGDPH2   | •     |
| Navigator Drill Positioning Handle 3              | SGDPH3   | •     |
| Navigator Drill Positioning Handle 4              | SGDPH4   | •     |
| Navigator Drill Positioning Handle 5              | SGDPH5   | •     |
| Certain Manual Bone Profiler, 3.4 mmD x 4.1 mmP   | IMMBP    | •     |
| Certain Manual Bone Profiler, 4.1 mmD x 5 mmP     | IMBP4    | •     |
| Certain Manual Bone Profiler, 5 mmD x 6 mmP       | IMBP5    | •     |
| Narrow Right Angle Large Hexed Driver Tip, 24 mmL | RASH3N   | •     |
| Open End Wrench                                   | CW100    | •     |
| Ratchet Wrench                                    | WR150    | •     |
| Handpiece Connector                               | MDR10    | •     |
| Posterior Large Hexed Driver, 17 mmL              | PHD02N   | •     |
| Ratchet Extension - Short, 6 mmL                  | RE100    | •     |
| Ratchet Extension - Long, 15 mmL                  | RE200    | •     |
| Navigator Certain Implant Mount, 3.4 mmD x 1(L)   | MSGIIM1  | •     |
| Navigator Certain Implant Mount, 3.4 mmD x 2(L)   | MSGIIM2  | •     |
| Navigator Certain Implant Mount, 3.4 mmD x 3(L)   | MSGIIM3  | •     |
| Navigator Certain Implant Mount, 3.4 mmD x 4(L)   | MSGIIM4  | •     |
| Navigator Certain Implant Mount, 4.1 mmD x 1(L)   | SGIIM41  | •     |
| Navigator Certain Implant Mount, 4.1 mmD x 2(L)   | SGIIM42  | •     |
| Navigator Certain Implant Mount, 4.1 mmD x 3(L)   | SGIIM43  | •     |
| Navigator Certain Implant Mount, 4.1 mmD x 4(L)   | SGIIM44  | •     |
| Navigator Certain Implant Mount, 5 mmD x 1(L)     | SGIIM51  | •     |
| Navigator Certain Implant Mount, 5 mmD x 2(L)     | SGIIM52  | •     |
| Navigator Certain Implant Mount, 5 mmD x 3(L)     | SGIIM53  | •     |
| Navigator Certain Implant Mount, 5 mmD x 4(L)     | SGIIM54  | •     |

The following items are not included in the Navigator Surgical Kit, but can be ordered separately, as needed, through Zimmer Biomet Dental Customer Service:

| Description                            | Item No. |
|----------------------------------------|----------|
| Navigator Twist Drill, 2 mmD x Y(L)    | SGD2Y    |
| Navigator Twist Drill, 2 mmD x Z(L)    | SGD2Z    |
| Navigator Twist Drill, 2.75 mmD x Y(L) | SGD275Y  |
| Navigator Twist Drill, 2.75 mmD x Z(L) | SGD275Z  |
| Navigator Twist Drill, 3 mmD x Y(L)    | SGD3Y    |
| Navigator Twist Drill, 3 mmD x Z(L)    | SGD3Z    |

| Description                            | Item No. |
|----------------------------------------|----------|
| Navigator Twist Drill, 3.25 mmD x Y(L) | SGD325Y  |
| Navigator Twist Drill, 3.25 mmD x Z(L) | SGD325Z  |
| Navigator Twist Drill, 3.85 mmD x Y(L) | SGD385Y  |
| Navigator Twist Drill, 3.85 mmD x Z(L) | SGD385Z  |
| Navigator Twist Drill, 4.25 mmD x Y(L) | SGD425Y  |
| Navigator Twist Drill, 4.25 mmD x Z(L) | SGD425Z  |

## Parallel Walled Navigator Laboratory Kit

SGLKIT

Parallel Walled Navigator Laboratory Kit



| Description                                    | Item No. | SGLKIT |
|------------------------------------------------|----------|--------|
| Navigator Laboratory Tray                      | SGLTRAY  | •      |
| Navigator Certain Analog Mount, 3.4 mmD x 1(L) | MSGIAM1  | •      |
| Navigator Certain Analog Mount, 3.4 mmD x 2(L) | MSGIAM2  | •      |
| Navigator Certain Analog Mount, 3.4 mmD x 3(L) | MSGIAM3  | •      |
| Navigator Certain Analog Mount, 3.4 mmD x 4(L) | MSGIAM4  | •      |
| Navigator Certain Analog Mount, 4.1 mmD x 1(L) | SGIAM41  | •      |
| Navigator Certain Analog Mount, 4.1 mmD x 2(L) | SGIAM42  | •      |
| Navigator Certain Analog Mount, 4.1 mmD x 3(L) | SGIAM43  | •      |
| Navigator Certain Analog Mount, 4.1 mmD x 4(L) | SGIAM44  | •      |
| Navigator Certain Analog Mount, 5 mmD x 1(L)   | SGIAM51  | •      |
| Navigator Certain Analog Mount, 5 mmD x 2(L)   | SGIAM52  | •      |
| Navigator Certain Analog Mount, 5 mmD x 3(L)   | SGIAM53  | •      |
| Navigator Certain Analog Mount, 5 mmD x 4(L)   | SGIAM54  | •      |

### Tapered Implant Surgical Kits

#### QNTSK40 Premium Kit



#### SLIMKT Slim Kit - Certain

#### SLIMKTE Slim Kit - Certain And External Hex



| Description                             | Item No. | Premium Kit | Upgrade Kit | Starter (Basic) Kit | Certain/External Hex | Certain |
|-----------------------------------------|----------|-------------|-------------|---------------------|----------------------|---------|
|                                         |          | QNTSK40     | QNTSK40U    | QNTSK20             | SLIMKTE              | SLIMKT  |
| Quad Shaping Drill, 3.25 mmD x 8.5 mmL  | QSD3285  | •           | •           |                     | •                    | •       |
| Quad Shaping Drill, 3.25 mmD x 10 mmL   | QSD3210  | •           | •           | •                   | •                    | •       |
| Quad Shaping Drill, 3.25 mmD x 11.5 mmL | QSD3211  | •           | •           | •                   | •                    | •       |
| Quad Shaping Drill, 3.25 mmD x 13 mmL   | QSD3213  | •           | •           | •                   | •                    | •       |
| Quad Shaping Drill, 3.25 mmD x 15 mmL   | QSD3215  | •           | •           |                     |                      |         |
| Quad Shaping Drill, 4 mmD x 8.5 mmL     | QSD485   | •           | •           |                     | •                    | •       |
| Quad Shaping Drill, 4 mmD x 10 mmL      | QSD410   | •           | •           | •                   | •                    | •       |
| Quad Shaping Drill, 4 mmD x 11.5 mmL    | QSD411   | •           | •           | •                   | •                    | •       |
| Quad Shaping Drill, 4 mmD x 13 mmL      | QSD413   | •           | •           | •                   | •                    | •       |
| Quad Shaping Drill, 4 mmD x 15 mmL      | QSD415   | •           | •           |                     |                      |         |
| Quad Shaping Drill, 5 mmD x 8.5 mmL     | QSD585   | •           | •           |                     |                      |         |
| Quad Shaping Drill, 5 mmD x 10 mmL      | QSD510   | •           | •           | •                   |                      |         |
| Quad Shaping Drill, 5 mmD x 11.5 mmL    | QSD511   | •           | •           | •                   |                      |         |
| Quad Shaping Drill, 5 mmD x 13 mmL      | QSD513   | •           | •           | •                   |                      |         |
| Quad Shaping Drill, 5 mmD x 15 mmL      | QSD515   | •           | •           |                     |                      |         |
| Quad Shaping Drill, 6 mmD x 8.5 mmL     | QSD685   | •           | •           |                     |                      |         |
| Quad Shaping Drill, 6 mmD x 10 mmL      | QSD610   | •           | •           |                     |                      |         |
| Quad Shaping Drill, 6 mmD x 11.5 mmL    | QSD611   | •           | •           |                     |                      |         |
| Quad Shaping Drill, 6 mmD x 13 mmL      | QSD613   | •           | •           |                     |                      |         |
| Quad Shaping Drill, 6 mmD x 15 mmL      | QSD615   | •           | •           |                     |                      |         |
| Dense Bone Tap, 3.25 mmD x 8.5 mmL      | NTAP3285 | •           | •           |                     | •                    | •       |
| Dense Bone Tap, 3.25 mmD x 10 mmL       | NTAP3210 | •           | •           | •                   | •                    | •       |
| Dense Bone Tap, 3.25 mmD x 11.5 mmL     | NTAP3211 | •           | •           | •                   | •                    | •       |
| Dense Bone Tap, 3.25 mmD x 13 mmL       | NTAP3213 | •           | •           | •                   | •                    | •       |
| Dense Bone Tap, 3.25 mmD x 15 mmL       | NTAP3215 | •           | •           |                     |                      |         |
| Dense Bone Tap, 4 mmD x 8.5 mmL         | NTAP485  | •           | •           |                     | •                    | •       |
| Dense Bone Tap, 4 mmD x 10 mmL          | NTAP410  | •           | •           | •                   | •                    | •       |
| Dense Bone Tap, 4 mmD x 11.5 mmL        | NTAP411  | •           | •           | •                   | •                    | •       |
| Dense Bone Tap, 4 mmD x 13 mmL          | NTAP413  | •           | •           | •                   | •                    | •       |
| Dense Bone Tap, 4 mmD x 15 mmL          | NTAP415  | •           | •           |                     |                      |         |

# Surgical Kits

## Site Preparation Technologies

### Tapered Implant Surgical Kits (Continued)

|                                                          |           | Premium Kit | Upgrade Kit | Starter (Basic) Kit | Certain/External Hex | Certain |
|----------------------------------------------------------|-----------|-------------|-------------|---------------------|----------------------|---------|
| Description                                              | Item No.  | QNTSK40     | QNTSK40U    | QNTSK20             | SLIMKTE              | SLIMKT  |
| Dense Bone Tap, 5 mmD x 8.5 mmL                          | NTAP585   | •           | •           |                     |                      |         |
| Dense Bone Tap, 5 mmD x 10 mmL                           | NTAP510   | •           | •           | •                   |                      |         |
| Dense Bone Tap, 5 mmD x 11.5 mmL                         | NTAP511   | •           | •           | •                   |                      |         |
| Dense Bone Tap, 5 mmD x 13 mmL                           | NTAP513   | •           | •           | •                   |                      |         |
| Dense Bone Tap, 5 mmD x 15 mmL                           | NTAP515   | •           | •           |                     |                      |         |
| Dense Bone Tap, 6 mmD x 8.5 mmL                          | NTAP685   | •           | •           |                     |                      |         |
| Dense Bone Tap, 6 mmD x 10 mmL                           | NTAP610   | •           | •           |                     |                      |         |
| Dense Bone Tap, 6 mmD x 11.5 mmL                         | NTAP611   | •           | •           |                     |                      |         |
| Dense Bone Tap, 6 mmD x 13 mmL                           | NTAP613   | •           | •           |                     |                      |         |
| Dense Bone Tap, 6 mmD x 15 mmL                           | NTAP615   | •           | •           |                     |                      |         |
| Tapered Implant Dense Bone Tap Kit Tray                  | NTAPKT    | •           | •           | •                   |                      |         |
| Tapered Implant Depth Indicator, 3.25 mmD x 8.5 mmL 3ea. | NTDI3285* | •           | •           |                     | •                    | •       |
| Tapered Implant Depth Indicator, 3.25 mmD x 10 mmL       | NTDI3210  | •           | •           | •                   | •                    | •       |
| Tapered Implant Depth Indicator, 3.25 mmD x 11.5 mmL     | NTDI3211  | •           | •           | •                   | •                    | •       |
| Tapered Implant Depth Indicator, 3.25 mmD x 13 mmL       | NTDI3213  | •           | •           | •                   | •                    | •       |
| Tapered Implant Depth Indicator, 3.25 mmD x 15 mmL       | NTDI3215  | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 4 mmD x 8.5 mmL 3ea.    | NTDI485*  | •           | •           |                     | •                    | •       |
| Tapered Implant Depth Indicator, 4 mmD x 10 mmL          | NTDI410   | •           | •           | •                   | •                    | •       |
| Tapered Implant Depth Indicator, 4 mmD x 11.5 mmL        | NTDI411   | •           | •           | •                   | •                    | •       |
| Tapered Implant Depth Indicator, 4 mmD x 13 mmL          | NTDI413   | •           | •           | •                   | •                    | •       |
| Tapered Implant Depth Indicator, 4 mmD x 15 mmL          | NTDI415   | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 5 mmD x 8.5 mmL 3ea.    | NTDI585   | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 5 mmD x 10 mmL          | NTDI510   | •           | •           | •                   |                      |         |
| Tapered Implant Depth Indicator, 5 mmD x 11.5 mmL        | NTDI511   | •           | •           | •                   |                      |         |
| Tapered Implant Depth Indicator, 5 mmD x 13 mmL          | NTDI513   | •           | •           | •                   |                      |         |
| Tapered Implant Depth Indicator, 5 mmD x 15 mmL          | NTDI515   | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 6 mmD x 8.5 mmL 3ea.    | NTDI685   | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 6 mmD x 10 mmL          | NTDI610   | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 6 mmD x 11.5 mmL        | NTDI611   | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 6 mmD x 13 mmL          | NTDI613   | •           | •           |                     |                      |         |
| Tapered Implant Depth Indicator, 6 mmD x 15 mmL          | NTDI615   | •           | •           |                     |                      |         |
| Certain Countersink Drill, 4 mmD                         | ICD100    | •           | •           | •                   | •                    | •       |
| QSD/NTDI Surgical Tray w/Plastic Insert                  | QNTSKT    | •           |             | •                   |                      |         |
| Open End Wrench                                          | CW100     | •           |             | •                   | •                    |         |
| Ratchet Wrench                                           | WR150     | •           |             | •                   |                      |         |
| High Torque Indicating Ratchet Wrench                    | H-TIRW**  |             |             |                     | •                    | •       |
| Ratchet Extension (Short), 6 mmL                         | RE100     | •           |             | •                   | •                    | •       |
| Ratchet Extension (Long), 15 mmL                         | RE200     | •           |             | •                   |                      |         |

\*One (1) each for SLIMKT and SLIMKTE.

\*\*H-TIRW is manufactured by Elos Med Tech and distributed by Zimmer Biomet Dental - Packaged Individually.

## Tapered Implant Surgical Kits (Continued)

| Description                                                 | Item No. | Premium Kit | Upgrade Kit | Starter (Basic) Kit | Certain/External Hex | Certain |
|-------------------------------------------------------------|----------|-------------|-------------|---------------------|----------------------|---------|
|                                                             |          | QNTSK40     | QNTSK40U    | QNTSK20             | SLIMKTE              | SLIMKT  |
| Implant Depth Gauge, 20 mmL                                 | DP020    | •           |             |                     |                      |         |
| Narrow Posterior Hexed Driver, 17 mmL                       | PHD00N   | •           |             | •                   | •                    |         |
| Posterior Large Hexed Driver, 17 mmL                        | PHD02N   | •           |             | •                   | •                    | •       |
| Narrow Right Angle Large Hexed Driver Tip, 24 mmL           | RASH3N   | •           |             | •                   |                      |         |
| Direction Indicator, 2 mmD x 10 mmL 3ea.                    | DI100    | •           |             |                     |                      |         |
| Handpiece Connector                                         | MDR10    | •           |             |                     | •                    | •       |
| Implant Mount, 3.4 mmD x 15 mmL                             | MMC15    | •           |             |                     | •                    |         |
| ACT Pointed Starter Drill                                   | ACTPSD   | •           |             |                     | •                    | •       |
| ACT Reusable Twist Drill, 2 mmD x 10 mmL                    | ACT2010  | •           |             |                     |                      |         |
| ACT Reusable Twist Drill, 2 mmD x 15 mmL                    | ACT2015  | •           |             |                     | •                    | •       |
| ACT Reusable Twist Drill, 2 mmD x 20 mmL                    | ACT2020  | •           |             | •                   |                      |         |
| ACT Drill Extension                                         | ACTDE    | •           |             |                     | •                    | •       |
| Internal Connection Universal Placement Driver Tip (Short)  | IIPDTUS  | •           |             | •                   | •                    | •       |
| Certain Implant Placement Driver Tip (Short), 4.1, 5, 6 mmD | IIPDTS   | •           |             | •                   | •                    | •       |
| Certain Implant Placement Driver Tip (Short), 3.4 mmD       | IMPDTS   | •           |             | •                   | •                    | •       |
| Internal Connection Universal Placement Driver Tip (Long)   | IIPDTUL  | •           |             | •                   |                      |         |
| Certain Implant Placement Driver Tip (Long), 4.1, 5, 6 mmD  | IIPDTL   | •           |             | •                   |                      |         |
| Certain Implant Placement Driver Tip (Long), 3.4 mmD        | IMPDTL   | •           |             |                     |                      |         |
| Certain Universal Ratchet Extension Driver (Short), 6 mmL   | IRE100U  | •           |             |                     |                      |         |
| Certain Universal Ratchet Extension Driver (Long), 15 mmL   | IRE200U  | •           |             | •                   | •                    | •       |
| Certain Driver Replacement O-Rings, 10-pack                 | IRORDR   | •           |             | •                   |                      |         |

# Surgical Kits

## Site Preparation Technologies

### Parallel Walled Surgical Kits

PSKT40

Premium Kit



Premium Kit

Certain Premium Kit

Starter (Basic) Kit

| Description                                  | Item No. | PSKT40 | PSKT35 | PSKT01 |
|----------------------------------------------|----------|--------|--------|--------|
| Plastic Surgical Tray Organizer              | PTT300I  | •      | •      | •      |
| Plastic Tray Organizer for PTT300I           | PTT100   | •      |        |        |
| ACT Reusable Twist Drill, 2 mmD x 10 mmL     | ACT2010  | •      | •      |        |
| ACT Reusable Twist Drill, 2 mmD x 15 mmL     | ACT2015  | •      | •      | •      |
| ACT Reusable Twist Drill, 2 mmD x 20 mmL     | ACT2020  | •      | •      |        |
| ACT Reusable Twist Drill, 2.75 mmD x 10 mmL  | ACT2710  | •      | •      |        |
| ACT Reusable Twist Drill, 2.75 mmD x 15 mmL  | ACT2715  | •      | •      | •      |
| ACT Reusable Twist Drill, 2.75 mmD x 20 mmL  | ACT2720  | •      |        |        |
| ACT Reusable Twist Drill, 3.25 mmD x 10 mmL  | ACT3210  | •      | •      |        |
| ACT Reusable Twist Drill, 3 mmD x 15 mmL     | ACT3015  | •      | •      | •      |
| ACT Reusable Twist Drill, 3.85 mmD x 20 mmL  | ACT3820  | •      | •      | •      |
| ACT Reusable Twist Drill, 4.25 mmD x 8.5 mmL | ACT4285  | •      | •      |        |
| ACT Reusable Twist Drill, 4.25 mmD x 13 mmL  | ACT4213  | •      | •      | •      |
| ACT Reusable Twist Drill, 4.25 mmD x 18 mmL  | ACT4218  | •      | •      |        |
| ACT Reusable Twist Drill, 4.85 mmD x 8.5 mmL | ACT4885  | •      | •      |        |
| ACT Reusable Twist Drill, 4.85 mmD x 13 mmL  | ACT4813  | •      | •      |        |
| ACT Reusable Twist Drill, 4.85 mmD x 18 mmL  | ACT4818  | •      | •      |        |
| ACT Reusable Twist Drill, 5.25 mmD x 8.5 mmL | ACT5285  | •      | •      |        |
| ACT Reusable Twist Drill, 5.25 mmD x 13 mmL  | ACT5213  | •      | •      |        |
| ACT Reusable Twist Drill, 5.25 mmD x 18 mmL  | ACT5218  | •      | •      |        |
| ACT Drill Extension                          | ACTDE    | •      | •      | •      |
| Countersink/Pilot Drill, 5 mmD               | CD500    | •      | •      | •      |

## Parallel Walled Surgical Kits (Continued)

| Description                                                 | Item No. | Premium Kit | Certain Premium Kit | Starter (Basic) Kit |
|-------------------------------------------------------------|----------|-------------|---------------------|---------------------|
|                                                             |          | PSKT40      | PSKT35              | PSKT01              |
| Certain Countersink Drill, 4 mmD                            | ICD100   | •           | •                   | •                   |
| Countersink Depth Indicator for ICD100                      | ICDI100  | •           | •                   | •                   |
| Countersink Drill, 4.1 mmD                                  | CD100    | •           |                     | •                   |
| Countersink/Pilot Drill, 6 mmD                              | CD600    | •           | •                   |                     |
| Pilot Drill, 2-3 mmD                                        | PD100    | •           | •                   | •                   |
| Implant Dense Bone Tap, 3.25 mmD x 18 mmL                   | MTAP2    | •           | •                   |                     |
| Implant Dense Bone Tap, 4 mmD x 13 mmL                      | TAP413   | •           | •                   |                     |
| Implant Dense Bone Tap, 5 mmD x 13 mmL                      | TAP53S   | •           | •                   |                     |
| Implant Dense Bone Tap, 6 mmD x 13 mmL                      | TAP63S   | •           |                     |                     |
| Osseotite 2 Dense Bone Tap, 5 mmD x 7-13 mmL                | XTAP53S  | •           | •                   |                     |
| Osseotite 2 Dense Bone Tap, 6 mmD x 7-13 mmL                | XTAP63S  | •           |                     |                     |
| Internal Connection Universal Placement Driver Tip (Short)  | IIPDTUS  | •           | •                   | •                   |
| Internal Connection Universal Placement Driver Tip (Long)   | IIPDTUL  | •           | •                   | •                   |
| Certain Implant Driver Tip (Short), 3.4 mmD                 | IMPDTS   | •           | •                   |                     |
| Certain Implant Driver Tip (Long), 3.4 mmD                  | IMPDTL   | •           | •                   | •                   |
| Certain Implant Placement Driver Tip (Short), 4.1, 5, 6 mmD | IIPDTS   | •           | •                   | •                   |
| Certain Implant Placement Driver Tip (Long), 4.1, 5, 6 mmD  | IIPDTL   | •           | •                   | •                   |
| Certain Universal Ratchet Extension Driver (Short), 6 mmL   | IRE100U  | •           | •                   | •                   |
| Certain Universal Ratchet Extension Driver (Long), 15 mmL   | IRE200U  | •           | •                   | •                   |
| Narrow Small Hexed Hand Driver, 17 mmL                      | PHD00N   | •           |                     | •                   |
| Narrow Large Hexed Hand Driver, 17 mmL                      | PHD02N   | •           | •                   | •                   |
| Narrow Large Hexed Driver, Tip 24 mmL                       | RASH3N   | •           | •                   | •                   |
| Driver Tip Replacement O-Rings, 10-pack                     | IRORDR   | •           | •                   |                     |
| Handpiece Connector                                         | MDR10    | •           | •                   | •                   |
| Open End Wrench                                             | CW100    | •           |                     | •                   |
| Ratchet Wrench                                              | WR150    | •           | •                   | •                   |
| Ratchet Extension (Short), 6 mmL                            | RE100    | •           |                     | •                   |
| Ratchet Extension (Long), 15 mmL                            | RE200    | •           | •                   | •                   |
| Implant Depth Gauge, 20 mmL                                 | DP020    | •           | •                   | •                   |
| Implant Mount, 3.4 mmD x 15 mmL                             | MMC15    | •           |                     | •                   |
| Direction Indicator, 2 mmD x 10 mmL 3ea.                    | DI100    | •           | •                   | •                   |
| Stainless Steel Round Drill, 2 mmD                          | RD100    | •           | •                   | •                   |

# Surgical Drills

## Site Preparation Technologies

### W & H Surgical Drilling Unit And Osstell ISQ Module\*

The surgical drill system is ideal for implantology with its automatic thread cutter function, automatic torque control, and the Osstell ISQ module for measuring implants stability.



| Description                                                                       | Item No. |
|-----------------------------------------------------------------------------------|----------|
| <b>Drill Systems</b>                                                              |          |
| SI-1023 control unit, 230 V                                                       | 30288000 |
| SI-1015 control unit, 120 V                                                       | 30289000 |
| <b>Motors</b>                                                                     |          |
| EM-19 LC electric motor for SI-1015, w/LED & light contacts (LC), 1.8 m cable     | 30281000 |
| EM-19 electric motor for SI-1015, without LED+ & light contacts (LC), 1.8 m cable | 30185000 |
| EM-19 LC electric motor for SI-1015, w/light contacts (LC), 3.5 m cable           | 30281001 |
| <b>Handpieces</b>                                                                 |          |
| WS-75 L, contra-angle 20:1, mini LED+ for EM-19 LC motor                          | 30033000 |
| WS-56 L, contra-angle 1:1, mini LED+ for EM-19 LC motor                           | 30040000 |
| WS-56 L G, contra-angle 1:1, self-generating mini LED+                            | 30039000 |
| WS-75 L G, contra-angle 20:1, self-generating mini LED+                           | 30032000 |
| <b>Foot Controls</b>                                                              |          |
| S-NW wireless foot control w/CAN dongle; compatible w/SI-1023, SI-1015            | 30264000 |
| S-NW wireless foot control w/SPI dongle; compatible w/ SI-915 & Piezomed SA-320   | 30264001 |
| S-N2 foot control w/cable, compatible w/SI-1015, SI-1023                          | 30285000 |
| <b>Accessories</b>                                                                |          |
| W&H Osstell ISQ module, for SI-1015 ,SI-1023                                      | 30210000 |
| W&H Osstell Measuring Probe                                                       | 07721100 |
| Wifi dongle CAN-Bus; SI-1015, SI-1023                                             | 07759700 |
| Wifi dongle SPI; SI-915 & Piezomed SA-320                                         | 07795800 |
| Universal Motor Holder                                                            | 07721800 |
| W&H Osstell Test Peg                                                              | 07849900 |
| SmartPeg Mount                                                                    | 07460300 |
| <b>SmartPegs</b>                                                                  |          |
| W&H SmartPeg Type 1 - 100353, Pack of 5                                           | 08009800 |
| W&H SmartPeg Type 4 - 100350, Pack of 5                                           | 08010120 |
| W&H SmartPeg Type 15 - 100386, Pack of 5                                          | 08009860 |
| W&H SmartPeg Type 26 - 100425, Pack of 5                                          | 08009970 |
| W&H SmartPeg Type 27 - 100431, Pack of 5                                          | 08009980 |
| W&H SmartPeg Type 30 - 100436, Pack of 5                                          | 08010010 |
| W&H SmartPeg Type 32 - 100440, Pack of 5                                          | 08010020 |
| W&H SmartPeg Type 33 - 100442, Pack of 5                                          | 08010030 |
| W&H SmartPeg Type 34 - 100444, Pack of 5                                          | 08010040 |
| W&H SmartPeg Type 35 - 100446, Pack of 5                                          | 08010050 |
| W&H SmartPeg Type 45 - 100472, Pack of 5                                          | 08010150 |
| W&H SmartPeg Type 60 - 100543, Pack of 5                                          | 08010620 |

\*Not available in all markets

## Advanced Cutting Technology (ACT) Twist Drills

Features include:

- Two-flute design for efficient cutting
- Alternating lines and bands enhance depth mark visibility and drill orientation in the osteotomy
- Matte finish to reduce glare
- Non-Irrigated Reusable



| Implant Lengths | Drill Lengths | 2 mmD   | 2.3 mmD | 2.75 mmD | 3 mmD   | 3.15 mmD | 3.25 mmD | 3.85 mmD | 4.25 mmD | 4.85 mmD | 5.25 mmD |
|-----------------|---------------|---------|---------|----------|---------|----------|----------|----------|----------|----------|----------|
| 7 – 8.5 mm      | 8.5 mm        | —       | —       | —        | —       | —        | —        | —        | ACT4285  | ACT4885  | ACT5285  |
| 7 – 10 mm       | 10 mm         | ACT2010 | —       | ACT2710  | ACT3010 | ACT3110  | ACT3210  | ACT3810  | —        | —        | —        |
| 7 – 13 mm       | 13 mm         | —       | —       | —        | —       | —        | —        | —        | ACT4213  | ACT4813  | ACT5213  |
| 7 – 15 mm       | 15 mm         | ACT2015 | —       | ACT2715  | ACT3015 | ACT3115  | ACT3215  | ACT3815  | —        | —        | —        |
| 7 – 18 mm       | 18 mm         | —       | —       | —        | —       | —        | —        | —        | ACT4218  | ACT4818  | ACT5218  |
| 7 – 20 mm       | 20 mm         | ACT2020 | —       | ACT2720  | ACT3020 | ACT3120  | ACT3220  | ACT3820  | —        | —        | —        |

## Quad Shaping Drills (QSDs) For Tapered Implants

Features include:

- Quad Shaping Drills use stepped depth marks for accuracy in creating the osteotomy
- QSDs have end cutting ability
- QSDs are color-coded for platform recognition



| Drill Lengths | 3.25 mmD | 4 mmD  | 5 mmD  | 6 mmD  |
|---------------|----------|--------|--------|--------|
| 8.5 mm        | QSD3285  | QSD485 | QSD585 | QSD685 |
| 10 mm         | QSD3210  | QSD410 | QSD510 | QSD610 |
| 11.5 mm       | QSD3211  | QSD411 | QSD511 | QSD611 |
| 13 mm         | QSD3213  | QSD413 | QSD513 | QSD613 |
| 15 mm         | QSD3215  | QSD415 | QSD515 | QSD615 |

## Miscellaneous Drills

| Description                        | Disposable | Reusable |
|------------------------------------|------------|----------|
| Round Drill, 2 mm                  | DR100      | RD100    |
| ACT Pointed Starter Drill          | —          | ACTPSD   |
| Pilot Drill, 2-3 mm                | DP100      | PD100    |
| Countersink Drill - Certain 4.1 mm | —          | CD100    |
| Countersink/Pilot Drill, 5 mm      | DC500      | CD500    |
| Countersink/Pilot Drill, 6 mm      | DC600      | CD600    |



# Surgical Components

## Site Preparation Technologies

### Dense Bone Taps For Parallel Walled Implants



| Implant Lengths | 3.25 mmD       | 4 mmD  | 5 mmD    | 6 mmD    |
|-----------------|----------------|--------|----------|----------|
| 8.5 – 18 mm     | MTAP1 (18 mmL) | —      | —        | —        |
| 8.5 – 18 mm     | MTAP2 (27 mmL) | —      | —        | —        |
| 7 – 10 mm       | —              | TAP410 | —        | —        |
| 7 – 13 mm       | —              | TAP413 | —        | —        |
| 7 – 20 mm       | —              | TAP420 | —        | —        |
| 6.5/7 – 8.5 mm  | —              | —      | XTAP58S  | XTAP68S  |
| 6.5/7 – 13 mm   | —              | —      | XTAP53S  | XTAP63S  |
| 6.5/7 – 18 mm   | —              | —      | XTAP518S | XTAP618S |

### Dense Bone Taps For Tapered Implants



| Implant Lengths | 3.25 mmD | 4 mmD   | 5 mmD   | 6 mmD   |
|-----------------|----------|---------|---------|---------|
| 8.5 mm          | NTAP3285 | NTAP485 | NTAP585 | NTAP685 |
| 10 mm           | NTAP3210 | NTAP410 | NTAP510 | NTAP610 |
| 11.5 mm         | NTAP3211 | NTAP411 | NTAP511 | NTAP611 |
| 13 mm           | NTAP3213 | NTAP413 | NTAP513 | NTAP613 |
| 15 mm           | NTAP3215 | NTAP415 | NTAP515 | NTAP615 |

NTAPK

Dense Bone Tap Kit

NOTE: Kit includes all of the Dense Bone Taps listed above.



## Tapered Implant Depth Indicators

Features include:

- Tapered Implant Depth Indicators (NTDIs) verify the position of the implant platform prior to placement
- NTDIs are color-coded for platform recognition




| Implant Lengths | 3.25 mmD | 4 mmD   | 5 mmD   | 6 mmD   |
|-----------------|----------|---------|---------|---------|
| 8.5 mm          | NTDI3285 | NTDI485 | NTDI585 | NTDI685 |
| 10 mm           | NTDI3210 | NTDI410 | NTDI510 | NTDI610 |
| 11.5 mm         | NTDI3211 | NTDI411 | NTDI511 | NTDI611 |
| 13 mm           | NTDI3213 | NTDI413 | NTDI513 | NTDI613 |
| 15 mm           | NTDI3215 | NTDI415 | NTDI515 | NTDI615 |

Also available in a kit:


**NTDIK** Depth/Direction Indicator Kit

NOTE: Kit includes all of the indicators listed to the left.



## Direction Indicators

| Description                                                      | Item No. |
|------------------------------------------------------------------|----------|
| 2 mm/3 mmD x 10 mmL                                              | DI100    |
| 2 mm/3 mmD x 15 mmL                                              | DI200    |
| 2.3 mm/3 mmD x 10 mmL, Anodized Blue                             | DI2310   |
| 2.3 mm/3 mmD x 15 mmL, Anodized Gold                             | DI2315   |
| Certain Direction Indicator Kit                                  | IDIKT    |
| Certain Direction Indicator for 4 mmD, 5 mmD, and 6 mmD Implants | IDI001   |
| Certain Direction Indicator for 3.25 mmD Implants                | IDIM01   |



# Surgical Components

## Site Preparation Technologies

### Handpiece Bone Profilers (Packaged With Guide Pin)

| Description                           | Item No. |
|---------------------------------------|----------|
| Bone Profiler, 3.4 mmD x 3.8 mmP      | BP340    |
| Bone Profiler, 3/4, 4.1 mmD x 5 mmP   | BP450    |
| Bone Profiler, 3/4, 4.1 mmD x 6 mmP   | BP460    |
| Bone Profiler, 3/4, 4.1 mmD x 7.5 mmP | BP475    |
| Bone Profiler, 4/5, 5 mmD x 5 mmP     | BP550    |
| Bone Profiler, 4/5, 5 mmD x 6 mmP     | BP560    |
| Bone Profiler, 4/5, 5 mmD x 7.5 mmP   | BP575    |
| Bone Profiler, 5/6, 6 mmD x 6 mmP     | BP660    |
| Bone Profiler, 5/6, 6 mmD x 7.5 mmP   | BP675    |



### Low Profile Abutment Bone Profilers

| Description                                 | Item No. |
|---------------------------------------------|----------|
| Low Profile Abutment Bone Profiler, 3.4 mmD | LPCBP3   |
| Low Profile Abutment Bone Profiler, 4.1 mmD | LPCBP4   |
| Low Profile Abutment Bone Profiler, 5 mmD   | LPCBP5   |



### Certain Manual Bone Profilers

| Description                                                | Item No. |
|------------------------------------------------------------|----------|
| Certain Manual Bone Profiler, 3.4 mmD x 4.1 mmP            | IMMBP    |
| Certain Manual Bone Profiler, 3.4 mmD x 5 mmP              | IMMBP35  |
| Certain Manual Bone Profiler, 4.1 mmD x 5 mmP              | IMBP4    |
| Certain Manual Bone Profiler, 5 mmD x 6 mmP                | IMBP5    |
| Certain Low Profile Abutment Manual Bone Profiler, 3.4 mmD | ILPCMBP3 |
| Certain Low Profile Abutment Manual Bone Profiler, 4.1 mmD | ILPCMBP4 |
| Certain Low Profile Abutment Manual Bone Profiler, 5 mmD   | ILPCMBP5 |



### Bone Profiler Guide Pins

| Description                                            | Item No. |
|--------------------------------------------------------|----------|
| Certain Bone Profiler Guide Pin, 3.25 mmD 3-pack       | IBPMGP   |
| Certain Bone Profiler Guide Pin, 4.1, 5 & 6 mmD 3-pack | IBPGP    |



### Bone Profiler Organizers And Kits

| Description                                                                                                    | Item No. |
|----------------------------------------------------------------------------------------------------------------|----------|
| Bone Profiler Organizer Box, 4.1 mmD                                                                           | BPKITB   |
| Bone Profiler Organizer Box, 5 & 6 mmD                                                                         | BPAKTB   |
| Bone Profiler Organizer Kit                                                                                    | BB001    |
| Bone Profiler Kit - For Use With 4.1 mmD Implants<br>Contains: BP450, BP460, BP475, BB001                      | BPKIT    |
| Bone Profiler Kit - For Use With 5 mmD or 6 mmD Implants<br>Contains: BP550, BP560, BP575, BP660, BP675, BB001 | BPAKT    |



## Tapered Osteotomes

Available in straight and angled designs for use with Tapered Implants

| Description                          |                | Item No. |
|--------------------------------------|----------------|----------|
| Straight Osseotite Tapered           | Osteotome Kit  | NTOST0   |
| Angled Osseotite Tapered             | Osteotome Kit  | NTOST0A  |
| Straight Osseotite Tapered, 3.25 mmD | Osteotome      | NTOST3   |
| Straight Osseotite Tapered, 4 mmD    | Osteotome      | NTOST4   |
| Straight Osseotite Tapered, 5 mmD    | Osteotome      | NTOST5   |
| Straight Osseotite Tapered, 6 mmD    | Osteotome      | NTOST6   |
| Angled Osseotite Tapered, 3.25 mmD   | Osteotome      | NTOST3A  |
| Angled Osseotite Tapered, 4 mmD      | Osteotome      | NTOST4A  |
| Angled Osseotite Tapered, 5 mmD      | Osteotome      | NTOST5A  |
| Angled Osseotite Tapered, 6 mmD      | Osteotome      | NTOST6A  |
| Straight Osseotite Tapered           | Osteotome Tray | NTOSTTRS |
| Angled Osseotite Tapered             | Osteotome Tray | NTOSTTRA |



## Summers Osteotomes

Available for use with Parallel Walled Implants

| Description                       | Item No. |
|-----------------------------------|----------|
| Summers Osteotome Kit #1-4        | OST00    |
| Summers Osteotome Kit #1-5        | OST10    |
| Summers Osteotome Kit #1-5, FS    | OST20    |
| Summers Osteotome #1              | OST01    |
| Summers Osteotome #2              | OST02    |
| Summers Osteotome #3              | OST03    |
| Summers Osteotome #4              | OST04    |
| Summers Osteotome #5              | OST05    |
| Summers Osteotome for Future Site | OSTFS    |
| Summers Osteotome System Tray     | OSTTR    |



## Sinus Elevators

| Description                          | Item No. |
|--------------------------------------|----------|
| Sinus Elevation Kit                  | SEKT1    |
| Titanium Nitride Coated Bone Carrier | TIBC1    |
| Small Membrane Elevator              | ME100    |
| Medium Membrane Elevator             | ME200    |
| Large Membrane Elevator              | ME300    |
| Large Titanium Plugger               | TIPL1    |
| Small Titanium Plugger               | TIPS1    |
| Sinus Elevation Tray                 | SEKTR    |



## Trephine Burs

| Description                                                | Item No. |
|------------------------------------------------------------|----------|
| Trephine Bur, 2 mm Inside Diameter/3 mm Outside Diameter   | TRE02    |
| Trephine Bur, 4 mm Inside Diameter/5 mm Outside Diameter   | TRE04    |
| Trephine Bur, 5 mm Inside Diameter/6 mm Outside Diameter   | TRE05    |
| Trephine Bur, 6 mm Inside Diameter/7 mm Outside Diameter   | TRE06    |
| Trephine Bur, 8 mm Inside Diameter/9.3 mm Outside Diameter | TRE08    |



D = Diameter • P = Restorative Platform

# Miscellaneous Instrumentation

## Site Preparation Technologies

| Description                              | Item No. |
|------------------------------------------|----------|
| Mallet Stainless Steel                   | MALL1    |
| Implant Depth Gauge                      | DP020    |
| Ratchet Extension, 6 mmL                 | IRE100U  |
| Ratchet Extension, 15 mmL                | IRE200U  |
| Titanium Elevator                        | TE003    |
| Manual Interchangeable Driver Tip Handle | MIDTH    |
| Tissue Measuring Post                    | TMP80    |



### H-TIRW High Torque Indicating Ratchet Wrench\*

Features include:

- The H-TIRW doubles as a conventional ratchet wrench to place implants and as a torque indicator to provided a visual assessment of torque
- Clearly indicates forward and reverse torque of 50-90 Ncm without the need for additional components
- Does not require lubrication or periodic re-calibration. Calibration may be checked by verifying that the arrow is positioned at the first scale mark (0 Ncm)
- Consisting of only two parts, the H-TIRW is easy to prepare for sterilization



\* L-TIRW and L-TIWRK are manufactured by Elos Medtech AB and distributed by Zimmer Biomet Dental.

| Description                                  |                  | Item No. |
|----------------------------------------------|------------------|----------|
| ACT Drill Extension                          |                  | ACTDE    |
| Certain Implant Placement Driver Tip - Long  |                  | IIPDTUL  |
| Certain Implant Placement Driver Tip - Short |                  | IIPDTUS  |
| Gelb Radiographic Depth Gauge Kit - 2 mmD    |                  | XDG00    |
| includes:                                    | 2 mmD x 13 mmL   | XDG13    |
|                                              | 2 mmD x 20 mmL   | XDG20    |
| Gelb Radiographic Depth Gauge Kit - 2.3 mmD  |                  | XDG01    |
| includes:                                    | 2.3 mmD x 13 mmL | XDG2313  |
|                                              | 2.3 mmD x 20 mmL | XDG2320  |
| Handpiece Connector                          |                  | MDR10    |
| Tissue Punch, 4.1 mm                         |                  | TE003    |
| Tissue Punch, 5 mm                           |                  | TP005    |
| Tissue Punch, 6 mm                           |                  | TP006    |



### BellaTek Encode Impression System

The BellaTek Encode Impression System eliminates the need for an implant-level impression while delivering a patient-specific restoration that has the appropriate margin height and natural emergence contours for the patient. This accuracy is achieved from an impression of the BellaTek Encode Healing Abutment. Codes embedded on the occlusal surface of the BellaTek Encode Healing Abutment communicate the collar height, implant hex-orientation, platform diameter and interface of the implant. Because the tissue is healed at the time of impression, the margin is designed to have the appropriate placement and contour. The result is a patient-specific, anatomically designed, final BellaTek Abutment.

### Efficient Workflow For The Entire Team

The process is simplified with the BellaTek Encode Healing Abutment!

#### Surgeon

- Efficient, streamlined interoffice processes simplify treatment for the referring dentist.
- Cutting-edge technologies create an improved and more simple treatment process, differentiating the practice to referring dentists and more importantly, to patients.
- The BellaTek Encode Impression System makes it easier for your referral base and may increase treatment acceptance.

#### Laboratory

- Potential new customers may lead to increased crown and bridge business.
- There is no need to create a cast, which results in less steps in the treatment process, reducing overhead.
- This unique branding opportunity may grow the volume of your business.

### The Result: A Highly Aesthetic BellaTek Definitive Abutment

#### Restorative Clinician

- Comfort: There is no need to use impression copings, resulting in a less invasive impression procedure for more comfort.
- Fewer Visits: The intraoral scan can be taken by the specialist at the surgical release visit, eliminating a restorative appointment and resulting in less visits to the dentist's office.
- Aesthetic Outcomes: Abutments designed specifically for the patient for better aesthetic outcomes.



BellaTek Encode Healing Abutment




BellaTek Abutments

Titanium  
Nitrate  
Coating


## Restorative Components

### Color-Coding

Certain Implants and all interfacing components are color-coded by restorative platform diameter for easy identification and selection.




Occlusal markings provide scanner information for abutment design.



### BellaTek Encode

| Restorative Platform | Emergence Profile | Collar Height |         |         |         |
|----------------------|-------------------|---------------|---------|---------|---------|
|                      |                   | 3 mmH         | 4 mmH   | 6 mmH   | 8 mmH   |
| 3.4 mmD              | 3.8 mmP           | IEHA343       | IEHA344 | IEHA346 | IEHA348 |
|                      | 5 mmP             | IEHA353       | IEHA354 | IEHA356 | IEHA358 |
| 4.1 mmD              | 4.1 mmP           | IEHA443       | IEHA444 | IEHA446 | IEHA448 |
|                      | 5 mmP             | IEHA453       | IEHA454 | IEHA456 | IEHA458 |
|                      | 6 mmP             | IEHA463       | IEHA464 | IEHA556 | IEHA468 |
| 5 mmD                | 5.6 mmP           | IEHA553       | IEHA554 | IEHA556 | IEHA558 |
|                      | 6 mmP             | IEHA563       | IEHA564 | IEHA566 | IEHA568 |
| 6 mmD                | 6.8 mmP           | IEHA663       | IEHA664 | IEHA666 | IEHA668 |



### EP One-Piece

| Restorative Platform | Emergence Profile | Collar Height |         |         |         |         |
|----------------------|-------------------|---------------|---------|---------|---------|---------|
|                      |                   | 2 mmH         | 3 mmH   | 4 mmH   | 6 mmH   | 8 mmH   |
| 3.4 mmD              | 3.4 mmP           | ISMHA32       | ISMHA33 | ISMHA34 | ISMHA36 | —       |
|                      | 3.8 mmP           | IMHA32        | IMHA33  | IMHA34  | IMHA36  | —       |
|                      | 5 mmP             | IMHA352       | IMHA353 | IMHA354 | IMHA356 | —       |
| 4.1 mmD              | 4.1 mmP           | ISHA42        | ISHA43  | ISHA44  | ISHA46  | ISHA48  |
|                      | 5 mmP             | ITHA52        | ITHA53  | ITHA54  | ITHA56  | ITHA58  |
|                      | 6 mmP             | —             | ITHA63  | ITHA64  | ITHA66  | ITHA68  |
| 5 mmD                | 5 mmP             | ISWHA52       | ISWHA53 | ISWHA54 | ISWHA56 | ISWHA58 |
|                      | 5.6 mmP           | IWTH52        | IWTH53  | IWTH54  | IWTH56  | IWTH58  |
|                      | 6 mmP             | IWTH562       | IWTH563 | IWTH564 | IWTH566 | IWTH568 |
| 6 mmD                | 6 mmP             | ISWHA62       | ISWHA63 | ISWHA64 | ISWHA66 | ISWHA68 |
|                      | 6.8 mmP           | IWTH62        | IWTH63  | IWTH64  | IWTH66  | IWTH68  |




### PreFormance Posts\*

**Material:** PEEK: Titanium Alloy Connection

**Indications:**

- Single- and multiple-unit cement-retained provisional restorations
- Minimum interarch space of 7 mm
- Maximum angulation of 15°
- Intraoral use limited to 180 days
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of occlusal forces in lateral and protrusive excursions and masticatory forces to minimize the load on the provisional restorations
- Guided soft-tissue healing of single and multiple-unit restorations of integrated implants



| Restorative Platform | Emergence Profile | Collar Height |                |          |                |
|----------------------|-------------------|---------------|----------------|----------|----------------|
|                      |                   | 4 mmH         |                | 6 mmH    |                |
|                      |                   | Straight      | 15° Pre-Angled | Straight | 15° Pre-Angled |
| 3.4 mmD              | 3.8 mmP           | IMPFP34       | IMPAPF34       | IMPFP36  | IMPAPF36       |
|                      | 5 mmP             | IMPFP354      | IMPAPF36       | IMPFP356 | IMPAPF356      |
| 4.1 mmD              | 5 mmP             | IPFP454       | IPAPF454       | IPFP456  | IPAPF456       |
| 5 mmD                | 5 mmP             | IWPFP554      | IPAPF554       | IWPFP556 | IPAPF556       |
| 6 mmD                | 6 mmP             | IWPFP664      | IPAPF664       | IWPFP666 | IPAPF666       |

### Temporary Cylinders\*

**Indications:**

- Single- and multiple-unit screw-retained provisional restorations
- Minimum interarch space of 6 mm
- Maximum angulation of 15°
- Intraoral use limited to 180 days for the PreFormance Temporary Cylinders
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of prosthesis support to minimize the load on the provisional abutments
- Occlusal loading of single- and multiple-unit restorations of integrated implants for guided tissue healing

### PreFormance Temporary Cylinders

**Material:** PEEK - Titanium Alloy Connection

| Restorative Platform | Description |           |
|----------------------|-------------|-----------|
|                      | Hexed       | Non-Hexed |
| 3.4 mmD              | IMPFTC34    | IMPFTC32  |
| 4.1 mmD              | IPFTC41     | IPFTC42   |
| 5 mmD                | IWPFTC51    | IWPFTC52  |
| 6 mmD                | IWPFTC61    | IWPFTC62  |



### Titanium Temporary Cylinders

**Material:** Titanium Alloy Connection

| Restorative Platform | Description |           |
|----------------------|-------------|-----------|
|                      | Hexed       | Non-Hexed |
| 3.4 mmD              | IMMTCS1     | —         |
| 4.1 mmD              | IITCS41     | IITCS42   |
| 5 mmD                | IWTCS51     | IWTCS52   |
| 6 mmD                | IWTCS61     | IWTCS62   |



\*All Provisional Abutments (PreFormance Posts and Temporary Cylinders) include the corresponding Titanium Hexed Screw.

D = Diameter • P = Emergence Profile • H = Collar Height


## Restorative Components

### GingiHue Posts

**Material:** Titanium Alloy

**Indications:**

- Single- and multiple-unit porcelain fused-to-metal restorations
- Areas of thin tissue where gingival discoloration is possible
- Preparable to follow gingival contours
- Conventional crown and bridge procedures
- Maximum angulation correction of 15°
- Minimum interarch space of 7 mm




| Restorative Platform | Emergence Profile | Collar Height |                |          |                |
|----------------------|-------------------|---------------|----------------|----------|----------------|
|                      |                   | 2 mmH         |                | 4 mmH    |                |
|                      |                   | Straight      | 15° Pre-Angled | Straight | 15° Pre-Angled |
| 3.4 mmD              | 3.8 mmP           | IMAP32G       | IMPAP32G       | IMAP34G  | IMPAP34G       |
|                      | 5 mmP             | IMAP352G      | IMPAP352G      | IMAP354G | IMPAP354G      |
| 4.1 mmD              | 5 mmP             | IAPP452G      | IPAP452G       | IAPP454G | IPAP454G       |
|                      | 6 mmP             | IAPP462G      | IPAP462G       | IAPP464G | IPAP464G       |
| 5 mmD                | 5 mmP             | IWPP552G      | IPAP552G       | IWPP554G | IPAP554G       |
|                      | 6 mmP             | IWPP562G      | IPAP562G       | IWPP564G | IPAP564G       |
| 6 mmD                | 6 mmP             | IWPP662G      | IPAP662G       | IWPP664G | IPAP664G       |

### UCLA Abutments

**Material:** Machined Gold Alloy Cylinder With Plastic Unitube

**Indications:**

- Cement- and screw-retained restorations
- Single- and multi-unit restorations
- Laboratory fabricated custom abutments
- Minimum interarch space of 6 mm
- For aesthetic restorations when tissue is limited
- Maximum angulation correction of 30°



| Description                                      | Material | 3.4 mmD   | 4.1 mmD   | 5 mmD     | 6 mmD     |
|--------------------------------------------------|----------|-----------|-----------|-----------|-----------|
| Hexed                                            | Gold     | IMUCG1C*  | IGUCA1C*  | IWGA51C*  | IWGA61C*  |
| Non-Hexed (With Gold-Tite® Large Diameter Screw) | Gold     | IMUCG2C   | IGUCA2C   | IWGA52C   | IWGA62C   |
| Non-Hexed (With Ti Large Diameter Screw)         | Gold     | IMUCG2T   | IGUCA2T   | IWGA52T   | IWGA62T   |
| Non-Hexed (With Gold-Tite Large Diameter Screw)  | Castable | IMUCC2C** | IUNAB2C** | IWPC52C** | IWPC62C** |
| Non-Hexed (With Ti Large Diameter Screw)         | Castable | MUCC2T**  | IUNAB2T** | IWPC52T** | IWPC62T** |

\* UCLA Gold Hexed Abutments do not include the screw. The final Gold-Tite Screw (IUNIHG) or the Titanium Screw (IUNIHT) must be ordered separately.

\*\* Not available in the US.

| Gold Alloy Cylinder Formulation |                                     |
|---------------------------------|-------------------------------------|
| Melting Range                   | 1400–1490°C (2550–2710°F)           |
| Solidus                         | 1400°C                              |
| Liquidus                        | 1490°C                              |
| CTE                             | 12.2 x 10 <sup>-6</sup> °K at 500°C |

| Gold Alloy Cylinder Melting Temperatures |                                     |                                     |
|------------------------------------------|-------------------------------------|-------------------------------------|
|                                          | Fahrenheit                          | Celsius                             |
| Melting Range                            | 2552-2732°F                         | 1400-1500°C                         |
| Solidus                                  | 2552-2579°F                         | 1400-1415°C                         |
| Liquidus                                 | 2714-2732°F                         | 1490-1500°C                         |
| CTE                                      | 12.2 x 10 <sup>-6</sup> °K at 932°F | 12.2 x 10 <sup>-6</sup> °K at 500°C |

# Definitive Restorations

## Restorative Components

### GingiHue Post And UCLA Abutment Components

#### Pick-Up Impression Coping\*

| Restorative Platform | Emergence Profile | Item No.                     |
|----------------------|-------------------|------------------------------|
| 3.4 mm               | 3.4 mm            | IMIC30                       |
|                      | 3.8 mm            | IMIC33                       |
|                      | 5 mm              | IMIC35                       |
| 4.1 mm               | 4.1 mm            | IIIC41<br>IIIC42 (non-hexed) |
|                      | 5 mm              | IIIC12                       |
|                      | 6 mm              | IIIC60                       |
|                      | 7.5 mm            | IIIC75                       |
| 5 mm                 | 5 mm              | IWIP50                       |
|                      | 5.6 mm            | IWIP55                       |
|                      | 6 mm              | IWIP56                       |
|                      | 7.5 mm            | IWIP57                       |
| 6 mm                 | 6 mm              | IWIP60                       |
|                      | 6.8 mm            | IWIP66                       |
|                      | 7.5 mm            | IWIP67                       |



#### Twist Lock™ Impression Coping\*\*

| Restorative Platform | Emergence Profile | Item No. |
|----------------------|-------------------|----------|
| 3.4 mm               | 3.8 mm            | IMIT33   |
|                      | 5 mm              | IMIT35   |
| 4.1 mm               | 4.1 mm            | IIIC44   |
|                      | 5 mm              | IIIC45   |
|                      | 6 mm              | IIIC46   |
|                      | 7.5 mm            | IIIC47   |
| 5 mm                 | 5 mm              | IWIT50   |
|                      | 5.6 mm            | IWIT55   |
|                      | 6 mm              | IWIT56   |
|                      | 7.5 mm            | IWIT57   |
| 6 mm                 | 6 mm              | IWIT60   |
|                      | 6.8 mm            | IWIT66   |
|                      | 7.5 mm            | IWIT67   |



#### Laboratory Analog

| Restorative Platform | Item No. |
|----------------------|----------|
| 3.4 mmD              | IMMILA   |
| 4.1 mmD              | IILA20   |
| 5 mmD                | IILAW5   |
| 6 mmD                | IILAW6   |



#### Additional Components

| Description                | Item No. |
|----------------------------|----------|
| Gold-Tite Screw            | IUNIHG   |
| Try-In Screws (5-pack)     | IUNITS   |
| Laboratory Abutment Holder | ILTAH57  |



\*Requires the use of Hex Driver [PHD02N](#), [PHD03N](#) or Driver Tip [RASH3N](#), [RASH8N](#)


\*\*Requires the use of Impression Coping Driver [ICD00](#)

## Low Profile Abutments

**Material:** Titanium Alloy

**Indications:**

- Single- and multiple-unit screw-retained restorations
- Minimum interarch space of 7.5 mm
- Minimum tissue height of 1 mm
- Angle correction up to 30°



| Restorative Platform | Collar Height | Straight (Non-Hexed) | Straight (Hexed) | 17° Pre-Angled (Hexed) | 30° Pre-Angled (Hexed) |
|----------------------|---------------|----------------------|------------------|------------------------|------------------------|
| 3.4 mmD              | 1 mm          | ILPC341U             | ILPC341          | —                      | —                      |
|                      | 2 mm          | ILPC342U             | ILPC342          | ILPAC3217              | —                      |
|                      | 3 mm          | ILPC343U             | ILPC343          | —                      | ILPAC3330              |
|                      | 4 mm          | ILPC344U             | ILPC344          | ILPAC3417              | —                      |
|                      | 5 mm          | —                    | —                | —                      | ILPAC3530              |
| 4.1 mmD              | 1 mm          | ILPC441U             | ILPC441          | —                      | —                      |
|                      | 2 mm          | ILPC442U             | ILPC442          | ILPAC4217              | —                      |
|                      | 3 mm          | ILPC443U             | ILPC443          | —                      | ILPAC4330              |
|                      | 4 mm          | ILPC444U             | ILPC444          | ILPAC4417              | —                      |
|                      | 5 mm          | —                    | —                | —                      | ILPAC4530              |
| 5 mmD                | 1 mm          | LPC541U              | ILPC541          | —                      | —                      |
|                      | 2 mm          | ILPC542U             | ILPC542          | ILPAC5217              | —                      |
|                      | 3 mm          | ILPC543U             | ILPC543          | —                      | ILPAC5330              |
|                      | 4 mm          | ILPC544U             | ILPC544          | ILPAC5417              | —                      |
|                      | 5 mm          | —                    | —                | —                      | ILPAC5530              |

## Low Profile Abutment Components

| Description                                                                                                                  | 4.8 mm Restorative Platform |
|------------------------------------------------------------------------------------------------------------------------------|-----------------------------|
| Pick-Up Impression Coping (hexed)         | LPCPIC1                     |
| Pick-Up Impression Coping (non-hexed)     | LPCPIC2                     |
| Twist Lock Impression Coping (hexed)      | LPCTIC1                     |
| Twist Lock Impression Coping (non-hexed)  | LPCTIC2                     |
| QuickBridge® Impression Coping                                                                                               | LPCRIC                      |
| Laboratory Analog                                                                                                            | LPCLA                       |
| Temporary Cylinder (hexed)                                                                                                   | LPCTC1                      |
| Temporary Cylinder (non-hexed)                                                                                               | LPCTC2                      |
| Healing Cap                               | LPCHC                       |
| PreFormance Temp. Cylinder (hexed)                                                                                           | LPCPTC1                     |
| PreFormance Temp. Cylinder (non-hexed)                                                                                       | LPCPTC2                     |
| QuickBridge                                                                                                                  | LPCQB                       |
| QuickBridge Replacement Cap                                                                                                  | LPCQBCAP                    |
| Distal Extension                                                                                                             | LPCDE                       |

| Description                                                                                           | 4.8 mm Restorative Platform |
|-------------------------------------------------------------------------------------------------------|-----------------------------|
| Gold Cylinder (hexed)                                                                                 | LPCGC1                      |
| Gold Cylinder (non-hexed)                                                                             | LPCGC2                      |
| Castable Cylinder (hexed)                                                                             | LPCCC1                      |
| Castable Cylinder (non-hexed)                                                                         | LPCCC2                      |
| Laser Welded Cylinder                                                                                 | LPCTUCA                     |
| Waxing Screw                                                                                          | LPCWS                       |
| Retaining Screw  | LPCGSH, LPCTSH              |
| Lapping Tool                                                                                          | LPCAMI                      |
| Polishing Protector                                                                                   | LPCPP                       |

| Abutment Compatibility |                |                    |
|------------------------|----------------|--------------------|
| Component              | Hexed Abutment | Non-Hexed Abutment |
| Hexed                  | Yes            | No                 |
| Non-Hexed              | Yes            | Yes                |

 = Requires the use of Hex Driver [PHD02N](#), [PHD03N](#) or Driver Tip [RASH3N](#), [RASH8N](#)  
 = Requires the use of Impression Coping Driver [ICD00](#)

D = Diameter • P = Emergence Profile • H = Collar Height

# Overdenture Restorations

## Restorative Components

### LOCATOR R-Tx Abutments

#### Material:

Abutment: Titanium Alloy With DuraTec® Titanium Carbon Nitride Coating

Housing: Titanium Alloy with pink anodization.

Inserts: Nylon

#### Indications:

- Tissue supported removable overdentures on 2-4 implants
- Partially edentulous overdentures with one or more implants
- Limited interarch distance
- Angle correction with up to 60° between divergent implants

| Collar Height | Restorative Platform |          |
|---------------|----------------------|----------|
|               | 3.4 mmP              | 4.1 mmP  |
| 1 mmH         | 30300-01             | 30301-01 |
| 2 mmH         | 30300-02             | 30301-02 |
| 3 mmH         | 30300-03             | 30301-03 |
| 4 mmH         | 30300-04             | 30301-04 |
| 5 mmH         | 30300-05             | 30301-05 |
| 6 mmH         | 30300-06             | 30301-06 |



### LOCATOR R-Tx Abutment Components

| Description                                    | Item No.  |
|------------------------------------------------|-----------|
| Zero Retention Insert - Grey, 4-pack           | 30001-01* |
| Low Retention Insert - Blue, 4-pack            | 30002-01* |
| Medium Retention Insert - Pink, 4-pack         | 30003-01* |
| High Retention Insert - Clear, 4-pack          | 30004-01* |
| Extra High Retention Insert - Green, 4-pack    | 30005-01  |
| Processing Insert, 4-pack                      | 30012-01  |
| Denture Attachment Processing Assembly, 4-pack | 30013-01* |
| 3.35 mm Abutment Analog, 4-pack                | 30014-01  |
| 4 mm Abutment Analog, 4-pack                   | 30015-01  |
| 5 mm Abutment Analog, 4-pack                   | 30016-01  |
| Impression Coping, 4-pack                      | 30017-01  |
| Processing Spacer, 4-pack                      | 30018-01  |
| Block-Out Spacer, 20-pack                      | 08514*    |



\* Included with each LOCATOR R-Tx Abutment.

### LOCATOR R-Tx Abutment Tools

| Description                    | Item No. |
|--------------------------------|----------|
| Retention Insertion Tool       | 30021-01 |
| Hexed Driver Short, 7 mmL      | PHD02N   |
| Hexed Driver Long, 24 mmL      | PHD03N   |
| Hexed Driver Tip Short, 24 mmL | RASH3N   |
| Hexed Driver Tip Long, 30 mmL  | RASH8N   |



LOCATOR R-Tx May Not Be Available In All Markets.

## LOCATOR Abutments

### Material:

Abutment: Titanium Alloy With Gold-Colored Titanium Nitride Coating

Housing: Titanium Alloy

Inserts: Nylon

### Indications:

- Tissue supported removable overdentures on 2-4 implants
- Partially edentulous overdentures with 1 or more implants
- Limited interarch distance
- Angle correction with up to 40° between divergent implants

| Collar Height | Restorative Platform |         |
|---------------|----------------------|---------|
|               | 3.4 mmP              | 4.1 mmP |
| 1 mmH         | IMLOA001             | ILOA001 |
| 2 mmH         | IMLOA002             | ILOA002 |
| 3 mmH         | IMLOA003             | ILOA003 |
| 4 mmH         | IMLOA004             | ILOA004 |
| 5 mmH         | IMLOA005             | ILOA005 |
| 6 mmH         | IMLOA006             | ILOA006 |



## LOCATOR Abutment Components

| Description                                               | Item No. |
|-----------------------------------------------------------|----------|
| Standard Replacement Insert - White, 4-pack               | LARMS*   |
| Extra Light Retention Insert - Blue, 4-pack               | LAELM*   |
| Light Retention Replacement Insert - Pink, 4-pack         | LLRMS*   |
| Extended Range Inserts - Green, 4-pack                    | LAERM    |
| Extra Light Retention Extended Range Insert - Red, 4-pack | LELARM   |
| Replacement Housing                                       | LOAH*    |
| Replacement Housing Kit                                   | LORHK    |
| Block-Out Spacer, 20-pack                                 | 08514*   |
| Impression Coping (one time use)                          | LAIC1    |
| Laboratory Analog                                         | LALA1    |



\* Included with each LOCATOR Abutment.

## LOCATOR Abutment Tools

| Description                                              | Item No.  |
|----------------------------------------------------------|-----------|
| Core Tool/Abutment Driver, 3-in-1 instrument             | LCTDR1    |
| Male Extraction Tool                                     | LOA8397** |
| Abutment Driver Tip Short, 24 mmL                        | LOADT4    |
| Abutment Driver Tip Large, 30 mmL                        | LOADT9    |
| 0.05 in/1.25 mm Hex Driver Tip 24 mmL - used with LCTDR1 | RASH4     |
| 0.05 in/1.25 mm Hex Driver Tip 30 mmL - used with LCTDR1 | RASH9     |

\*\*Part of the Core Tool/Abutment Driver LCTDR1.

D = Diameter • P = Restorative Platform • L = Length • H = Collar Height

## Restorative Components

## Abutment Screws



|               | Gold-Tite Screw                                                                       | Titanium Screw                                                                        | Large Gold-Tite Screw                                                      | Large Titanium Screw |
|---------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------------------------------------------------------------------------|----------------------|
| Screws        | IUNIHG                                                                                | IUNIHT                                                                                | ILRGHG                                                                     | ILRGHT               |
| Drivers       | PHD02N, PHD03N                                                                        | PHD02N, PHD03N                                                                        | PHD02N, PHD03N                                                             | PHD02N, PHD03N       |
| Driver Tips   | RASH3N, RASH8N                                                                        | RASH3N, RASH8N                                                                        | RASH3N, RASH8N                                                             | RASH3N, RASH8N       |
| Torque        | 20 Ncm                                                                                | 20 Ncm                                                                                | 20 Ncm                                                                     | 20 Ncm               |
| For Use With: | UCLA Abutment,<br>GingiHue Post,<br>GingiHue Post 15°,<br>Hexed Temporary<br>Cylinder | UCLA Abutment,<br>GingiHue Post,<br>GingiHue Post 15°,<br>Hexed Temporary<br>Cylinder | Non-Hexed UCLA Abutment,<br>Non-Hexed Temporary Cylinder,<br>BellaTek Bars |                      |

## Retaining Screws



|               | Low Profile Abutment<br>Gold-Tite Screw | Low Profile Abutment<br>Titanium Screw |
|---------------|-----------------------------------------|----------------------------------------|
| Screws        | LPCGSH                                  | LPCTSH                                 |
| Drivers       | PHD02N, PHD03N                          | PHD02N, PHD03N                         |
| Driver Tips   | RASH3N, RASH8N                          | RASH3N, RASH8N                         |
| Torque        | 10 Ncm                                  | 10 Ncm                                 |
| For Use With: | Low Profile Abutment                    | Low Profile Abutment                   |


Waxing Screws  
And Try-In Screws

|               | Laboratory Waxing Screw                                                      | Hexed Try-In Screw (5-pack)                                 |
|---------------|------------------------------------------------------------------------------|-------------------------------------------------------------|
| Screws        | IWSU30                                                                       | IUNITS                                                      |
| Drivers       | PHD02N, PHD03N                                                               | PHD02N, PHD03N                                              |
| Driver Tips   | RASH3N, RASH8N                                                               | RASH3N, RASH8N                                              |
| For Use With: | Pick-Up and<br>Impression Copings,<br>Implant Level<br>Laboratory Procedures | Hexed UCLA Abutment,<br>GingiHue Post,<br>GingiHue Post 15° |

## Restorative Components

BellaTek Bars and Frameworks are created with the technology to achieve a truly passive fit. This technology does not require any wax or resin design from the technician. Because these are milled from one solid piece of lightweight titanium, the potential for weakness caused by soldering or laser welded joints is eliminated. You receive a completed superstructure ready for prosthesis processing. No polishing or retention creation necessary. Zimmer Biomet Dental can even tap holes for attachments. Copymilled Bars are also available.

### Hybrid Bar Design #1




- Fixed prosthesis
- Metal lingual
- Low buccal finish line
- Titanium tissue surface

### Hybrid Bar Design #2



- Fixed prosthesis
- Metal lingual
- Low buccal finish line
- Titanium tissue surface

### Wraparound Bar



- Fixed prosthesis
- Acrylic resin wraparound design
- Adjustments possible to accommodate for changing tissue contours
- Easy to repair and modify the acrylic resin portion of the prosthesis

### Freeform Bar



- Fixed prosthesis
- Acrylic resin wraparound design

### DOLDER® Egg Shape Bar



- Removable overdenture design
- Anterior clips allow for anterior/posterior rotation of the overdenture
- Provides lateral stability to the overdenture

### DOLDER U Shape Bar



- Removable overdenture design
- Anterior clips allow for anterior/posterior rotation of the overdenture
- Provides lateral stability to the overdenture

### Hader Bar



- Removable overdenture design with clips that engage undercuts
- Provides lateral stability to the overdenture
- 4 mm between implants is needed for clips

### Primary Bar




- Most secure removable overdenture design
- May be used with attachments or secondary casting
- Bar can be designed with a 2° – 6° taper

### Combination Primary Bar With Hader Or DOLDER Designs



- Combines a milled bar with Hader or DOLDER Designs
- Provides excellent lateral stability to the overdenture
- Bar can be designed with a 2° – 6° taper

## Design Features



### Copymilled Cobalt Chrome Framework



- Fixed screw-retained prosthesis, allowing veneering of user friendly CoCr porcelains with CTE  $14.1 \times 10^{-6} K^{-1}$
- High biocompatible Cobalt Chrome Alloy, certified for orthopedic surgical implants

### Copymilled Titanium Framework



- Ideal for fixed screw-retained prosthesis. Scanned from a lab designed resin pattern
- Available in Titanium Alloy for acrylic Finishing or Commercially Pure titanium for Titanium porcelain veneer

## Restorative Components

NCATD0

Contra Angle Torque Driver Kit\*



NPSDK0

Prosthetic Driver Kit\*



NCATD0C

Contra Angle Torque Driver Kit Certain\*



L-TIRWK

Low Torque Indicating Ratchet Wrench Kit\*\*



| Description                                   | Item No.  | Contra Angle Torque Kit | Contra Angle Torque Kit Certain | Prosthetic Driver Kit | Low Torque Indicating Kit |
|-----------------------------------------------|-----------|-------------------------|---------------------------------|-----------------------|---------------------------|
|                                               |           | NCATD0                  | NCATD0C                         | NPSDK0                | L-TIRWK                   |
| Narrow Posterior Hex Driver, 17 mmL           | PHD02N    | •                       | •                               | •                     |                           |
| Narrow Posterior Square Driver, 17 mmL        | PSQD0N    | •                       |                                 | •                     |                           |
| Contra Angle Torque Driver Body               | HTD-C     | •                       | •                               | •                     |                           |
| Narrow Right Angle Hex Driver Tip, 24 mmL     | RASH3N    | •                       | •                               | •                     | •                         |
| Narrow Right Angle Hex Driver Tip, 30 mmL     | RASH8N    | •                       | •                               | •                     |                           |
| Right Angle Abutment Driver Tip - Steel       | RASA3     | •                       | •                               | •                     |                           |
| Narrow Right Angle Square Driver Tip, 24 mmL  | RASQ3N    | •                       |                                 | •                     |                           |
| Narrow Right Angle Square Driver Tip, 30 mmL  | RASQ8N    | •                       |                                 | •                     |                           |
| Contra Angle Torque Controller, 10 Ncm*       | TD-10     |                         |                                 | •                     |                           |
| Contra Angle Torque Controller, 20 Ncm*       | TD-20     | •                       | •                               | •                     |                           |
| Contra Angle Torque Controller, 35 Ncm*       | TD-35     | •                       |                                 | •                     |                           |
| Contra Pre-action Handle                      | CATDH     | •                       | •                               | •                     |                           |
| Prosthetic System Driver Tray                 | PSDT1     | •                       | •                               | •                     |                           |
| Narrow Impression Coping Driver               | ICD00     |                         |                                 | •                     |                           |
| Posterior Abutment Driver, 17 mmL             | PAD00     |                         |                                 | •                     |                           |
| Standard Abutment Driver, 24 mmL              | PAD24     |                         |                                 | •                     |                           |
| Posterior Slotted Screw Driver, 17 mmL        | PSD00     |                         |                                 | •                     |                           |
| Standard Slotted Screw Driver, 24 mmL         | PSD01     |                         |                                 | •                     |                           |
| Narrow Standard Slotted Square Driver, 24 mmL | PSQD1N    |                         |                                 | •                     |                           |
| Narrow Standard Hex Driver, 24 mmL            | PHD03N    |                         |                                 | •                     |                           |
| Right Angle Slotted Driver Tip, 30 mmL        | RASD1     |                         |                                 | •                     |                           |
| Low Torque Indicating Ratchet Wrench          | L-TIRW**  |                         |                                 |                       | •                         |
| Restorative Torque Tray                       | RTI2035TR |                         |                                 |                       | •                         |

\* NCATD0, NCATD0C and NPSDK0 are manufactured by NSK America Corp/NSK dental LLC and distributed by Zimmer Biomet Dental.

\*\*L-TIRW and L-TIWRK are manufactured by Elos Medtech AB and distributed by Zimmer Biomet Dental.

The cordless Prosthodontic Screwdriver has torque control, which enables precise tightening of retaining screws. This cordless screwdriver reduces the need for traditional mechanical rotary ratchet wrench systems.



### Features Include:

- Torque control of 8-40 Ncm in 1 Ncm steps for precise tightening of retaining screws
- An automatic stop when the preset torque is reached
- Cordless operation with constant drive speed
- A small contra-angle head for better visibility of the treatment site
- A contra-angle that can be rotated 360°
- A clockwise and counter clockwise operation
- A sterilizable contra-angle

### Prosthodontic Screwdriver IA-400

#### Includes:

- A handpiece drive IA-40H with keypad and rechargeable Li-Ion battery (not sterilizable)
- A contra-angle IA-80 with a hexagon chucking system
- A charging station with country specific cables available

| Description                         | Item No. |
|-------------------------------------|----------|
| IA-400 Prosthodontic Screw Driver   | 16934000 |
| IA-400 Charging Station             | 05338301 |
| IA-80 Contra-Angle Hand Piece, 80:1 | 16338001 |
| IA-40H Handpiece Drive              | 16834000 |
| Operation Lever                     | 05937000 |
| Power Cord US, CAN, J               | 05456100 |
| Power Cord UK, IRL                  | 05457700 |
| Power Cord EU                       | 05464400 |
| Power Cord AUS, NZ                  | 05467600 |
| Rechargeable Battery, Li-Ion        | 05451300 |

# Drivers, Driver Tips And Miscellaneous Restorative Products

38

## Restorative Components

|               |  |  |  |  |  |  |  |
|---------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|               | Abutment                                                                          |                                                                                   | LOCATOR                                                                           | Large Hex                                                                         |                                                                                     | Abutment                                                                            | LOCATOR                                                                             |
| Length        | Driver                                                                            | Tip                                                                               | Tip only                                                                          | Driver                                                                            | Tip                                                                                 | Driver Only                                                                         | Tip Only                                                                            |
| 17 mm         | PAD00                                                                             | —                                                                                 | —                                                                                 | PHD02N                                                                            | —                                                                                   | ICD00                                                                               | —                                                                                   |
| 24 mm         | PAD24                                                                             | RASA3                                                                             | LOADT4                                                                            | PHD03N                                                                            | RASH3N                                                                              | —                                                                                   | RASH4                                                                               |
| 30 mm         | —                                                                                 | —                                                                                 | LOADT9                                                                            | —                                                                                 | RASH8N                                                                              | —                                                                                   | RASH9                                                                               |
| For Use With: | Low Profile Abutments                                                             |                                                                                   | LOCATOR Abutments                                                                 | Healing Abutments, Abutment Screws, Certain Cover Screws and Retaining Screws     |                                                                                     | Twist Lock Transfer LOCATOR Abutments with Impression Copings                       | LOCATOR Core Tools                                                                  |

| Description                                                        | Item No. |
|--------------------------------------------------------------------|----------|
| Lapping Tool For Castable UCLA Abutment With Guide, 3.4 mmD        | IMLT150  |
| Lapping Tool For Castable UCLA Abutment With Guide, 4.1 mmD, 5 mmD | ILT150   |
| Reamer and Handle                                                  | IRH600   |
| Certain QuickSeat® Activator Tool                                  | IQSA01   |
| Polishing Protector, UCLA Abutment, 3.4 mmD                        | IPPM1    |
| Polishing Protector, UCLA Abutment, 4.1 mmD                        | IPPIA3   |
| Polishing Protector, UCLA Abutment, 5 mmD                          | IPPIA5   |
| Polishing Protector, UCLA Abutment, 6 mmD                          | IPPIA6   |
| Polishing Protector, UCLA Abutment, 3.4 mmD, 4.1 mmD, 5mmD, 6 mmD  | ILTAH57  |
| Low Torque Indicating Ratchet Wrench*                              | L-TIRW   |
| Low Torque Indicating Ratchet Wrench Kit*                          | L-TIRWK  |
| Standard ISO1797 Adapter                                           | C9980    |
| Short ISO1797 Adapter                                              | C9981    |
| Prosthetic Angle Guide Kit, 15°, 25°, 35°                          | AG900    |
| Manual Interchangeable Driver Tip Handle                           | MIDTH    |
| Tissue Measuring Post                                              | TMP80    |



\* L-TIRW and L-TIWRK are manufactured by Elos Medtech AB and distributed by Zimmer Biomet Dental.

D = Diameter

For comprehensive information, please see Regenerative Product Brochure or visit [www.zimmerbiometdental.com](http://www.zimmerbiometdental.com).





**ZIMMER BIOMET**  
Your progress. Our promise.

Contact us at 1-800-342-5454 or visit

[zimmerbiometdental.com](http://zimmerbiometdental.com)

Zimmer Biomet Dental  
Global Headquarters  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Tel: +1-561-776-6700  
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. LOCATOR, LOCATOR R-Tx, R-Tx, Color and Shape of Retention Inserts, and DuraTec are registered trademarks of ZEST IP Holdings, LLC. Implantmed is manufactured by W&H Dentalwerk Bürmoos GmbH. Implantmed and W&H are registered trademarks of W&H Dentalwerk Bürmoos GmbH. Dolder is a registered trademark of Dr. Charlotte Fleischer-Dolder and Dr. Fritz Dolder. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions.

This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental.

ZBINST1260 REV D 04/20 ©2020 Zimmer Biomet. All rights reserved.

